



**University of
Zurich^{UZH}**

**Zurich Open Repository and
Archive**

University of Zurich
University Library
Strickhofstrasse 39
CH-8057 Zurich
www.zora.uzh.ch

Year: 2013

Generating whole bacterial genome sequences of low-abundance species from complex samples with IMS-MDA

Seth-Smith, H M B ; Harris, S R ; Scott, P ; Parmar, S ; Marsh, P ; Unemo, M ; Clarke, I N ; Parkhill, J ; Thomson, N R

Abstract: The study of bacterial populations using whole genome sequencing is of considerable scientific and clinical interest. However, obtaining bacterial genomic information is not always trivial: the target bacteria may be difficult-to-culture or uncultured, and may be found within samples containing complex mixtures of other contaminating microbes and/or host cells, from which it is very difficult to derive robust sequencing data. Here we describe our procedure to generate sufficient target genomic DNA for whole bacterial genome sequencing, from clinical samples of the difficult-to-culture, obligate intracellular pathogen *Chlamydia trachomatis*, without the need for culture. Our protocol combines Immunomagnetic Separation (IMS) for targeted bacterial enrichment with Multiple Displacement Amplification (MDA) for whole genome amplification, followed by high throughput sequencing. Compared to other techniques which might be used to generate such data, IMS-MDA is an inexpensive, low-technology and highly transferable process, which provides amplified genomic DNA for sequencing from target bacteria in under 5 hours, with little hands-on time.

DOI: <https://doi.org/10.1038/nprot.2013.147>

Posted at the Zurich Open Repository and Archive, University of Zurich

ZORA URL: <https://doi.org/10.5167/uzh-86263>

Journal Article

Accepted Version

Originally published at:

Seth-Smith, H M B; Harris, S R; Scott, P; Parmar, S; Marsh, P; Unemo, M; Clarke, I N; Parkhill, J; Thomson, N R (2013). Generating whole bacterial genome sequences of low-abundance species from complex samples with IMS-MDA. *Nature Protocols*, 8(12):2404-2412.

DOI: <https://doi.org/10.1038/nprot.2013.147>

1 Generating whole bacterial genome sequences of low- 2 abundance species from complex samples with IMS-MDA

3
4
5 Helena M.B. Seth-Smith^{1,2*}, Simon R. Harris¹, Paul Scott¹, Surendra Parmar³, Peter Marsh⁴,
6 Magnus Unemo⁵, Ian N. Clarke⁶, Julian Parkhill¹ and Nicholas R. Thomson^{1*}
7

8 ¹ Pathogen Genomics, The Wellcome Trust Sanger Institute, Wellcome Trust Genome Campus, Hinxton,
9 Cambridgeshire, UK

10 ² Current address: Functional Genomics Center Zürich, Universität Zürich, Winterthurerstrasse 190, 8057-
11 Zürich, Switzerland

12 ³ Clinical Microbiology and Public Health Laboratory, Health Protection Agency, Addenbrooke's Hospital,
13 Cambridge, UK

14 ⁴ Public Health England Laboratory Southampton, Southampton General Hospital, Southampton, UK

15 ⁵ WHO Collaborating Centre for Gonorrhoea and Other STIs, National Reference Laboratory for Pathogenic
16 Neisseria, Örebro University Hospital, Örebro, Sweden

17 ⁶ Molecular Microbiology Group, Faculty of Medicine, University of Southampton, Southampton General
18 Hospital, Southampton, UK
19

20	Helena M.B. Seth-Smith*	hss@fgcz.uzh.ch
21	Simon R. Harris	sh16@sanger.ac.uk
22	Paul Scott	ps6@sanger.ac.uk
23	Surendra Parmar	surendra.parmar@addenbrookes.nhs.uk
24	Peter Marsh	peter.marsh@uhs.nhs.uk
25	Magnus Unemo	magnus.unemo@orebroll.se
26	Ian N. Clarke	inc@soton.ac.uk
27	Julian Parkhill	parkhill@sanger.ac.uk
28	Nicholas R. Thomson*	nrt@sanger.ac.uk
29		

30 *Corresponding authors

31 Helena M.B. Seth-Smith: T +41 44 6358579; F +41 44 635 39 22

32 Nicholas R. Thomson: T +44 1223 494740; F +44 1223 494919
33

34 Keywords: Immunomagnetic Separation, Multiple Displacement Amplification, Whole
35 Genome Amplification, Fastidious, Difficult-to-culture, *Chlamydia trachomatis*, high-
36 throughput, antibodies, magnetic beads
37

1 Abstract/summary

2 The study of bacterial populations using whole genome sequencing is of considerable scientific
3 and clinical interest. However, obtaining bacterial genomic information is not always trivial: the
4 target bacteria may be difficult-to-culture or uncultured, and may be found within samples
5 containing complex mixtures of other contaminating microbes and/or host cells, from which it is
6 very difficult to derive robust sequencing data. Here we describe our procedure to generate
7 sufficient target genomic DNA for whole bacterial genome sequencing, from clinical samples of
8 the difficult-to-culture, obligate intracellular pathogen *Chlamydia trachomatis*, without the need
9 for culture. Our protocol combines Immunomagnetic Separation (IMS) for targeted bacterial
10 enrichment with Multiple Displacement Amplification (MDA) for whole genome amplification,
11 followed by high throughput sequencing. Compared to other techniques which might be used to
12 generate such data, IMS-MDA is an inexpensive, low-technology and highly transferable process,
13 which provides amplified genomic DNA for sequencing from target bacteria in under 5 hours,
14 with little hands-on time.

15

1 Introduction

2 Driven by plummeting sequencing costs, whole genome sequencing has revolutionised the
3 way we have been able to understand bacterial biology, pathogenesis, epidemiology,
4 genetics and evolution. Sequencing technologies have been applied to track the temporal
5 and geographic distribution of pathogens, and offer the precision to understand the true
6 nature of infection within a single person¹⁻⁷. Much of our understanding of the nature and
7 diversity of bacteria derives from species which can be easily, rapidly and selectively grown
8 under laboratory conditions. It is a challenge to generate genomic DNA of sufficient quantity
9 and quality for whole genome sequencing from uncultured, fastidious or difficult-to-culture
10 bacteria. Discarded clinical samples represent an important resource for studying bacteria,
11 although such samples may be complex, often containing low levels of the species of
12 interest among a multitude of contaminating bacteria and host cells.

13 We have developed an approach to enrich for specific bacteria and amplify sufficient
14 genomic DNA for whole genome sequencing, directly from complex or non-viable samples
15 without the need for culture⁸. Our methodology combines Immunomagnetic Separation
16 (IMS) with Multiple Displacement Amplification (MDA) and allows for access to the genomes
17 of specific bacteria even where the target species is present at low levels and with other
18 contaminating microbiota and host cells. To validate our methodology we amplified and
19 sequenced complete genomes of *C. trachomatis* strains present in discarded clinical
20 samples, representing the first time that whole bacterial genome sequences have been
21 generated directly from uncultured clinical samples⁸.

22 *Chlamydia trachomatis* is a pathogen of global importance, causing more than 100 million
23 cases of sexually transmitted chlamydial infection annually⁹ as well as the blinding ocular
24 disease trachoma¹⁰. There are several subspecies-level typing schemes for *C. trachomatis*,
25 which group strains according to variation in the highly-variable gene *ompA*, a panel of
26 housekeeping genes (MLST)¹¹⁻¹³, or fast-changing repetitive loci (variable number tandem
27 repeats, VNTRs)¹⁴. However, these methods lack the resolution required for detailed strain
28 tracking in the case of *ompA*-genotyping and MLST, and can be vulnerable to stochastic
29 change in the case of VNTRs. Moreover, recently it has been shown that recombination in *C.*
30 *trachomatis*, leading to genetic exchange and diversification unlinked to phylogeny, severely
31 limits the interpretation of data from these typing methodologies¹⁵⁻¹⁷. For *C. trachomatis*,
32 as for many other bacterial species, whole genome sequencing is the only technology that
33 can provide the resolution required to determine true relationships, in addition to the
34 accuracy and specificity required to differentiate closely related isolates typical of
35 monomorphic species⁴.

36 *C. trachomatis* is an ideal model bacterium for testing non-culture based sequencing
37 protocols: as an obligate intracellular pathogen it requires tissue culture for *in vitro* growth,
38 a method which is technically challenging, expensive and time consuming; some strains may

be recalcitrant to culture; and discarded clinical samples are available for research. While deep-sequencing of a clinical sample has been shown to be able to give some information about an infecting *C. trachomatis* strain¹⁸, this method is expensive and is not scalable or able to provide sufficient strain resolution. Several other technologies exist or are being developed for the depletion of host cells in clinical samples (including MolYsis Basic by Molzym and GeneRead Bacterial DNA Kit by Qiagen [not yet released]), particularly blood samples (discussed in ¹⁹), but these are of less use when other microbiota are present in the sample, masking data from the bacterium of interest. Also of note are DNA target enrichment technologies, where the sequence of a reference strain can be used to develop baits covering the whole genome, allowing the specific hybridisation of the target strain's DNA. Such methods have been shown to be effective in several studies on bacteriophage, bacteria and viruses ¹⁹⁻²³ and can be used on samples in lysis buffer, but they are highly expensive and time-consuming, requiring the initial design and purchase of a custom array of baits. Additionally, if the species under study has a variable genome, not all the variation in the genome will be captured using these methods.

The sequencing of genomes of single bacteria has been shown in several studies, isolated by micromanipulation, microfluidics or flow cytometry (reviewed in ^{24,25}). For some applications such as comparison of sequence variation between individual cells of the same species within a complex mix, methodologies such as these may offer some advantages, depending on the specific research goal. However, these techniques are highly time and resource consuming, with the isolation stages so far not targeted towards a particular bacterial species. Moreover it is very difficult to produce full coverage of the genomes of interest, which is an aim of our protocol, such that detailed genomic epidemiology can be performed. Further recent techniques include a "mini-metagenome" approach following automated cell sorting, and growth of individual bacteria from communities within gel microdroplets prior to amplification ^{26,27}, from which better genome coverage is achieved as a larger number of bacteria are subject to amplification.

The need for methods to target the genomes of difficult-to-culture bacteria was further demonstrated by a very recent paper²⁸, in which IMS and MDA were also used in combination to sequence *C. trachomatis* genomes from clinical samples. This method is largely equivalent to the one detailed here, with similar antibodies employed, although alternative magnetic beads were used, a DNase step was included to reduce host DNA contamination (a step which we found did not increase the success rate), and an additional DNA extraction step was used prior to MDA.

Development of the protocol

IMS is an established technique for enriching target bacteria from complex mixtures using antibodies attached to magnetic beads. In the past, IMS has been used with varied samples,

in order to concentrate the target bacteria and remove inhibitors from the sample, enabling clearer and more accurate nucleic acid amplification-based diagnosis²⁹⁻³². We chose to test IMS for bacterial enrichment prior to genome sequencing.

Our experience of clinical swab samples for *C. trachomatis* diagnosis has shown that the vast majority of samples do not carry sufficient quantities of target bacteria and their DNA to allow genome sequencing from the sample either directly, or following amplification⁸. Consequently we chose to follow IMS with MDA which amplifies high molecular weight DNA using ø29 polymerase and random hexamer primers^{33,34}. This combined IMS-MDA protocol can provide sufficient high quality genomic DNA for sequencing using high throughput technologies⁸. IMS-MDA is a rapid, flexible, low-technology, low cost (for consumable reagents, excluding quantification and sequencing costs) protocol with high potential for use with a multitude of sample types and bacterial species.

Overview of IMS-MDA

The IMS-MDA protocol is shown in **Figure 1** and follows the following stages: 1) Primary and magnetic bead-conjugated secondary antibodies are mixed. 2) After incubation, a magnet is used to retain the bound primary and secondary antibodies, and excess unbound primary antibody is removed. 3) The bound antibodies are mixed with the clinical sample. 4) After incubation, samples are washed to remove contaminating microbes and cells which are not bound to the antibodies. 5) Enriched target cells are retained by a magnet. 6) MDA is performed directly on the enriched bacterial sample. Subsequent to IMS-MDA, the following steps are used to generate genome sequence data: 7) The amplified DNA is quantified, and samples with sufficient target DNA are sequenced using a high throughput sequencing platform. 8) Bioinformatic analysis is used to determine the quality of the sequence data and the nature of the targeted bacterial genome sequence.

Applications of the method

Our work focussed on obtaining complete genome sequences of the difficult-to-culture, sexually transmitted bacterium *C. trachomatis*⁸. The samples we first used to validate this approach were cultured and purified *C. trachomatis*, after which we tested the protocol on discarded urogenital swab samples that had returned a positive diagnosis for *C. trachomatis* by a routine diagnostic nucleic acid amplification test (NAAT)^{8,35}. We have also successfully applied this technique to non-viable archived *C. trachomatis* samples, providing access to historical genome data that would have otherwise been lost, as the samples contained insufficient DNA even for sequencing at high coverage levels, without prior enrichment.

1 We believe that the IMS-MDA protocol would be equally applicable to any micro-organism
2 for which a suitable and specific primary antibody or aptamer exists or can be developed,
3 and where suitable primary specimen material is available. The sample may be viable or
4 non-viable, as long as the bacterium remains intact such that the antigen remains associated
5 with the genomic DNA. The low cost and high transportability of this Protocol means that
6 samples can be prepared even in resource poor laboratories, and the critical first stages of
7 IMS (stages 1-5 in **Figure 1**; 1-15 in the Procedure) could even be performed in the field.
8 IMS-MDA has advantages over target enrichment techniques in situations when a reference
9 genome sequence is not available, or when dealing with a highly variable genome, as *de*
10 *novo* sequence assembly is possible using the sequence data generated by IMS-MDA.

11

12 **Experimental design**

13 **Sample choice.** This technique is appropriate for samples which maintain intact bacteria.
14 Our study used clinical samples placed directly into *Chlamydia* Transport Medium (CTM).
15 Appropriate ethical permission may be required prior to work on clinical samples. In order
16 to comply with the Human Tissue Act (2004) as it applies in England, Wales and Northern
17 Ireland, human cells must be lysed on receipt and sequence data from patient DNA must be
18 discarded prior to sequence analysis (see details in Materials).

19 Clinical samples which involve swabs or urine samples placed directly into lysis buffer are
20 not suitable for this approach. With the general move of molecular diagnostic laboratories
21 towards the use of commercial NAATs which supply bespoke collection vessels containing
22 lysis/stabilisation buffer, such as those for Gen-Probe Aptima or Abbott m2000rt, future
23 studies will generally require prospective additional samples to be taken and collected in
24 suitable transport media as described above. Other possible samples could include any
25 complex mixture with sufficient intact target material. Samples can be concentrated by
26 centrifugation from larger volumes prior to resuspension in a suitable medium to keep the
27 cells intact (e.g. transport medium or PBS), or homogenised if required. This method has so
28 far been unsuccessful in generating complete genome sequences from urine samples,
29 perhaps due to the condition of the bacteria. The minimum number of bacteria required for
30 a successful outcome was difficult to precisely establish using *C. trachomatis*, which is not
31 easily quantifiable, but a concentration dependent effect was observed during the
32 validation of the procedure⁸. The affinity of the chosen antibody and the amount of
33 contaminating material may also affect the required minimum input.

34 **Choice of antibody.** Selection of an appropriate antibody is key for IMS-MDA, as it is the
35 critical component for the specific enrichment of the target organism. Bacterial surface
36 antigens such as lipopolysaccharides (LPS) and outer membrane proteins are common
37 antibody targets. Previously developed commercial diagnostic antibodies are excellent
38 reagents for IMS, as they are likely to have been extensively tested for cross reactivity. For

our study we used the commercially available anti-*Chlamydia* mouse IgG primary antibody (see Materials) which targets chlamydial LPS, present at approx. 34,000 molecules per bacterium³⁶, and binds to all serovars of *C. trachomatis*. This reagent proved ideal since it is inexpensive for the quantities used, was originally developed for diagnostic purposes and thus has been tested for cross-reactivity against many other microbial species possibly present in the urogenital tract including *Lactobacillus lactis*, *Mycoplasma* spp., *Neisseria gonorrhoeae* and *Gardnerella vaginalis* (IMAGEN *Chlamydia* booklet^{37,38}). We also tested a range of other anti-*Chlamydia* antibodies which provided comparable results: mouse monoclonal IgG2b anti-*C. trachomatis* LPS, mouse monoclonal IgG2a anti-*C. trachomatis* MOMP and rabbit polyclonal IgG anti-*C. trachomatis* (see Materials). Similarly, small peptide aptamers that bind to specific target molecules could be used for this purpose.

Magnetic bead selection. The primary antibodies can be attached to magnetic beads through use of specific secondary, magnetic bead-conjugated antibodies, including anti-mouse or anti-rabbit IgG antibodies. The primary antibody used to validate this protocol is not directly linked to a magnetic bead and so we used a secondary anti-mouse IgG sheep antibody. Alternatively, primary antibodies can be directly conjugated to activated magnetic beads (e.g. Dynabeads® MyOne™ Carboxylic Acid, Invitrogen).

Whole genome amplification. MDA was chosen as the most appropriate method for whole genome amplification (WGA). Since the development of this protocol, a new method of WGA has been developed: multiple annealing and looping-based amplification cycles (MALBAC), which reduces the observed amplification bias of the isothermal MDA through cycles of primer extension³⁹. While the resulting genome coverage may be more even, MALBAC currently produces more errors in nucleotide incorporation through the use of DNA polymerases with lower fidelity than the ø29 polymerase. Therefore we continue to recommend the use of MDA.

Bacterial genome quantitation by qPCR. A specific qPCR system is required to enable quantification of the target species' genomic DNA obtained by IMS-MDA. Our initial experiments used a SYBR Green system, but greater accuracy, sensitivity and reproducibility was obtained with the Taqman® system described in Jalal *et al.* for *C. trachomatis*³⁵, targeting the single copy chromosomal *ompA* gene. The information derived from qPCR directs the decision on which samples are suitable for downstream sequencing. Quantification of the total DNA resulting from the MDA reaction, which will include amplified DNA from contaminating organisms, is unnecessary and uninformative for our purpose. For *C. trachomatis*, we found that a minimum of 1,500,000 genome copies per µl was required to provide complete genome sequences with confidence in single nucleotide polymorphisms (SNP) calling.

Controls. Fresh buffer can be used as a negative control for contamination and be processed in parallel with experimental samples, although we had never found contamination to occur. Buffer containing the target organism can be used as a positive control, if available.

1 Positive and negative controls for MDA are described with the manufacturer's instructions.
2 The final result is determined as being positive or negative after analysis of the sequencing
3 data.

4 **Genome sequencing.** We used Illumina GAII and HiSeq platforms in our study, but the DNA
5 produced should be equally amenable to other sequencing technologies. For SNP analysis,
6 high-throughput technologies such as Illumina or Ion Torrent are recommended because of
7 their high and accurate sequence yield per sequencing run. We have used high levels of
8 multiplexing (up to 96) with IMS-MDA samples with no deleterious effects.

9 **Data analysis.** We mapped the sequence data to a known reference genome to determine
10 the extent, depth and evenness of coverage. This was performed using SMALT
11 (<http://www.sanger.ac.uk/resources/software/smalt/>), although similar programs are
12 available including bwa and SOAP ^{40,41}. Using our experimental design, bases and SNPs can
13 be called with accuracy when the mean depth of coverage is greater than 35x, as long as the
14 coverage coefficient of variation (CV = standard deviation/mean) is not greater than 1. MDA
15 can produce uneven coverage, and with CV values >0.5 manual checking of base calling is
16 recommended. The *C. trachomatis* genome is extremely stable in terms of gene content and
17 overall genome architecture ^{16,42-45}, meaning that mapping-based approaches allowed us to
18 accurately reconstruct more than 99% of the genome of samples meeting these coverage
19 criteria. For other bacteria, where the genome may contain accessory regions that would
20 not be accessible by mapping-based methods, *de novo* assembly would be required. In such
21 cases it is important to note that many assembly algorithms make the assumption of
22 relatively constant read-coverage levels across the sequenced genome, and may produce
23 poor results when used with sequence reads from samples that have been subject to IMS-
24 MDA. We successfully assembled large contigs from *C. trachomatis* IMS-MDA data using
25 SPAdes ⁴⁶, an assembly program designed for use with sequence reads from MDA-amplified
26 single-cell-derived DNA. For *de novo* assembly of IMS-MDA data we would therefore
27 recommend use of an algorithm designed for sequence data from amplified DNA, such as
28 velvet-sc ⁴⁷, IDBA-UD ⁴⁸ or SPAdes ⁴⁶. These programs also account for chimeric DNA
29 rearrangements that may occur during MDA, and which again may confound analyses with
30 assembly programs not designed for this specific purpose. In order to remove assembled
31 regions derived from contaminating DNA, we aligned the resulting contigs against a
32 reference *C. trachomatis* genome sequence using abacas ⁴⁹. This process is enabled by the
33 extremely conserved nature of the *C. trachomatis* genome and allows the identification and
34 resolution of inverted repeats, which may otherwise result in artefacts during assembly. The
35 resulting *C. trachomatis* genomes in our study were finally assembled into 2 – 21 contigs ⁸.
36 We were also able to identify one mixed infection within the samples tested, and to
37 separate out the sequences of both strains based on the relative coverage of the two
38 genomes ⁸.

39

1 **Limitations of the method**

2 IMS-MDA allows access to the genomic information in otherwise inaccessible samples and
3 our studies have shown that we can generate accurate whole genome sequences from 15-
4 30% of discarded clinical swab samples that have tested positive by a routine diagnostic
5 NAAT ⁸. The reasons for lack of success in some samples is unclear, but was unrelated to
6 choice of antibody or incubation conditions. Therefore it is assumed that the load of the
7 target bacterium and the integrity of the sample are keys to the success of the technique.
8 Clearly samples in which the genomic DNA is dissociated from the targeted antigen will not
9 provide sufficient DNA of interest to sequence. Careful selection of samples is essential to
10 allow the greatest likelihood of success.

1 Materials

2 REAGENTS

3 NaH₂PO₄×H₂O (Sigma-Aldrich, cat. no. S9638)

4 Na₂HPO₄×2H₂O (Sigma-Aldrich, cat. no. S3264)

5 NaCl (Sigma-Aldrich, cat. no. S3014)

6 Tween 20® (Sigma-Aldrich, cat. no. P9416)

7 Primary antibodies: IMAGEN *Chlamydia* (Oxoid, ThermoFisher, cat. no. K610111-2), mouse
8 monoclonal anti-*C. trachomatis* LPS IgG2b (MyBioSource, clone no. M4020310), mouse
9 monoclonal anti-*C. trachomatis* MOMP IgG2a (AbCam, Cambridge, UK, clone no. BIOD166)
10 and rabbit polyclonal anti-*C. trachomatis* IgG (MyBioSource, San Diego, California, USA, cat.
11 no. MBS221885)

12 Dynabeads® M-280 Sheep Anti-Mouse IgG (Invitrogen, cat. no. 11201D) or other
13 appropriate secondary antibody-coupled or activated Dynabeads®.

14 illustra GenomiPhi V2 DNA Amplification Kit (GE Healthcare, cat. no. 25-6600-30) or RepliG
15 (Qiagen, cat. no. 150023)

16 Clinical samples. We obtained discarded urogenital clinical samples in Remel M4RT
17 transport medium (ThermoFisher). **!CAUTION** Appropriate ethical permission might be
18 required for the use of discarded clinical samples. Our study was approved by the
19 appropriate National Research Ethics Service Committee. Additionally, we automatically
20 mapped all sequence reads to the human genome directly after the sequence generation,
21 and discarded human sequence data without further analysis. **! CAUTION:** Clinical samples
22 will very likely contain viable bacteria, therefore samples should be handled under
23 appropriate containment conditions (for *C. trachomatis* under containment level 2,
24 preferably in a class 2 biological safety cabinet) until the bacteria have been heat-killed.

25 TaqMan® Fast Advanced Master Mix (Applied Biosystems, Invitrogen, cat. no. 4444556) or
26 alternative qPCR reagents.

27 qPCR primers and probes. Ideally the qPCR target should be a single-copy locus on the
28 bacterial chromosome. The largest chromosome should be targeted if there is more than
29 one in the bacterium of interest as MDA amplifies smaller circular molecules such as
30 plasmids to higher copy numbers than larger molecules. If no genomic information on the
31 bacterium exists, it could be attempted to target specific 16S rRNA gene sequences. For *C.*
32 *trachomatis* quantification, we used the Taqman primer set targeting the chromosomal
33 single copy *ompA* gene F primer HJ-MOMP-1: 5' GACTTTGTTTCGACCGTGTT, R primer HJ-

1 MOMP-2: 5' ACARAATACATCAAARCGATCCCA, probe MOMP: 5' VIC-
2 ATGTTTACVAAAYGCGCTT³⁵ (Sigma Aldrich custom oligos).

3

4

5 EQUIPMENT

6 DynaMag™-2 Magnet (Invitrogen, cat. no. 12321D).

7 Shaking incubator (Innova 42, New Brunswick Scientific, Eppendorf, cat. no. M1335-0002
8 or equivalent) or platform rocker (STR6, Stuart Scientific, ScienceLab.com, cat. no. 65-286-
9 674 or equivalent). A tube rotator can also be used (SB2 with SB3/ 1, Stuart Scientific, cat.
10 no. SB2 and SB3/1 or equivalent).

11 Microcentrifuge (Eppendorf 5418, cat. no. FA-45-18-11 or equivalent).

12 2 ml Safe-lock tubes™ (Eppendorf, cat. no. 0030 120.094).

13 Vortex (Vortex-Genie 2, Scientific Industries, cat. no. SI-0266 or equivalent).

14 PCR tubes or plates and sealer, 0.2 ml (Thin-walled dome-capped strips, Thermo Scientific,
15 cat. no. AB -0451, Eppendorf Twin.tec, Eppendorf, cat. no. 0030133374 and Microseal A
16 film, BioRad, cat no. MSA5002 or equivalent).

17 PCR machine (e.g. MJ Research Tetrad 2 DNA Engineer cyclers, Bio-Rad, cat. no. PTC-0240G).

18 Real-time quantitative PCR machine (StepOne™ Real-Time PCR System, Life Technologies,
19 Invitrogen, cat. no. 4376600, or equivalent).

20 MicroAmp Fast Optical 96 well qPCR plate (Applied Biosystems, cat. no. 4311971).

21 MicroAmp Fast Optical adhesive film (Applied Biosystems, cat. no. 4346906).

22 Where required: class 2 biological safety cabinet for handling clinical samples or other
23 samples which may contain live containment level 2 organisms.

24

25 REAGENT SETUP

26 Isotonic PBS. Mix 0.16 g/l NaH₂PO₄·H₂O, 0.98 g/l Na₂HPO₄·2H₂O, 8.10 g/l NaCl; autoclave
27 and aliquot into sterile 100 ml bottles.

28 PBST (PBS-Tween) wash buffer. Add Tween 20® to isotonic PBS to 0.05% vol/vol. Prepare
29 fresh before use. **CRITICAL:** Lack of Tween 20 in the wash buffer can cause clumping of the
30 beads and lower the efficacy of the procedure.

1 qPCR standards. Serially diluted standards should be used alongside experimental samples
2 for approximating the number of target organism genome copies. For this purpose, either
3 use a known quantity of pure genomic DNA or, where necessary, generate, purify and
4 quantify a PCR product covering the locus targeted by the qPCR assay. For the *C.*
5 *trachomatis ompA* Taqman assay we generated standard PCR products using the primers f:
6 5'-CGGAATTGTGCATTTACGTG3'; r: 5'-CTACGCTGAGGACGGTAAGC3'.

7

Helena 30.8.13 12:33

Kommentar [1]: No, this is a one off, making and diluting the standards. I would consider it part of the setup.

1 Procedure

2 IMS of bacteria from complex samples – TIMING minimum 2.5h

3 **CRITICAL** This protocol can be used to process between 1 and 16 samples simultaneously.
4 For each sample to be processed, use 2 µl of Dynabeads®. Beads for all samples are
5 prepared in one pool in steps 1-9, prior to aliquotting them to the individual samples.

6 **1** Vortex the stock of Dynabeads®, remove the required volume and place in a 2 ml
7 tube.

8 **2** Add 0.5 ml PBST and resuspend the beads with a 1 sec vortex or by flicking the tube.
9 Place the tube in the DynaMag™-2 Magnet and leave for 2 min. Remove the liquid
10 while the tube is still held in the magnet, retaining the Dynabeads® in the tube.

11 ? TROUBLESHOOTING

12 **3** Remove the tube from the magnet and repeat step 2.

13 **4** Remove the tube from the magnet and add 0.5 ml PBST. Per sample to be processed,
14 add primary antibody to the equivalent of approx. 10^{12} IgG molecules (e.g. 0.25 µl of
15 IMAGEN *Chlamydia* per sample).

16 **5** Incubate the beads at 20°C shaking at 200 rpm for at least 1 h. Note that we observe
17 no decrease in performance using temperatures in the range of 4-30°C, incubation
18 time of up to 24 h, or with alternative mixing methods including rotation and
19 rocking.

20 **6** During the incubation, prepare the clinical samples by defrosting them if necessary.
21 **!CAUTION** Samples containing live bacteria must be handled under appropriate
22 containment conditions (see Materials).

23 **7** Perform a pulse spin in the microcentrifuge on the tube from Step 5.

24 **8** Place the tube in the magnet and leave for 2 min. Remove the liquid while the tube
25 is still held in the magnet, retaining the Dynabeads® in the tube.

26 **9** To remove any unbound primary antibody from the solution, perform step 2 twice
27 and resuspend in 50 µl PBST per sample to be processed.

28 **10** Place aliquots of the clinical samples (10-200 µl depending on the available volume)
29 in 2 ml tubes. Add 50 µl of the antibody-bound Dynabeads® from Step 9 to each
30 sample and mix by flicking the tubes.

31 ? TROUBLESHOOTING

32 **11** Perform step 5.

33 **12** Perform a pulse spin on the tubes.

34 **13** Place the tubes in the magnet and leave for 2 min. Remove the liquid while the tubes
35 are still held in the magnet, retaining the Dynabeads® in the tubes.

36 ? TROUBLESHOOTING

37 **14** To remove contaminating material from the target bacteria, perform Step 2 twice.

15 After the final wash remove as much of the wash buffer as possible, leaving just the Dynabeads® in the tubes.

PAUSE POINT Material can be stored at -20°C for at least 14 days prior to amplification. Material can be transported at this stage. To inactivate biological material for transport, samples can be heated to 95°C for 5 min and subsequently stored at -20°C for several months.

MDA – TIMING 2.5 h

16 Perform MDA using the illustra GenomiPhi V2 DNA Amplification Kit. Resuspend each of the Dynabeads® pellets from Step 15, which carry the template DNA, in 9 µl sample buffer, and transfer the full volume to a 0.2 ml well of a PCR plate. The Dynabeads® can remain in the samples throughout the reaction.

17 Complete the MDA according to the manufacturer's instructions, using the 95°C denaturation step to release DNA from the bacteria.

DNA Quantification and Genome Sequencing – TIMING several days up to weeks, depending on method used

18 Quantify the amount of target bacterium DNA produced by performing Taqman qPCR. Make sufficient master mix for the required number of samples and standards using TaqMan® Fast Advanced Master Mix in a reaction volume of 20 µl, with each primer at a final concentration of 300nM and probe at final concentration of 150nM.

19 Add 1 µl of each sample from Step 17 to 19 µl of master mix, in a 96 well qPCR plate and seal with adhesive film.

20 Perform the qPCR using standard cycling conditions shown in the table.

Cycle number	Uracil-N glycosylase incubation	Polymerase activation	Denature	Anneal / Extend
1	50 °C, 2 min	95 °C, 20 s		
2–41			95 °C, 1 s	60 °C, 20s

21 Analyse the results in the qPCR software, using the cycle threshold (C_t) values from the standards to calculate the amount of target DNA in the IMS-MDA samples.

? TROUBLESHOOTING

CRITICAL STEP For obtaining complete, accurate bacterial genome sequences, >1,500,000 genome copies per µl are required (data from *C. trachomatis*). Minimum

Helena 30.8.13 12:59

Kommentar [2]: This is really going to depend on the software used. In the StepOne software, the values are automatically calculated for standards and samples, and given for each sample, which is why I haven't elaborated.

Zlotoryński, Eytan 30.8.13 10:41

Kommentar [3]: Au, this is the Step to explain (see also comment 3): Please explain the parameter to be recorded in the genomic/dilution standards qPCR reactions, and instruct to prepare a function from all the standards against which the experimental samples can be quantified.

requirement for other bacteria and samples should be empirically determined. It is not recommended to proceed to sequencing with samples with <500,000 genome copies per μ l.

22 Use the remaining material from Step 17 (19 μ l) for high throughput sequencing, up to 96-plex. (Illumina GAII or HiSeq, or Ion Torrent technology is suitable to generate a high yield of accurate sequence data.)

23 Analyse the resulting sequence reads by mapping to a known reference genome using appropriate software to determine the depth and evenness of coverage, such as SMALT, bwa or SOAP. At least 99% of the genome should have reads to map against it, with the mean depth of coverage greater than 35 \times , although these values should be verified for each organism and sequencing method used. Calculate the coverage coefficient of variation (CV = standard deviation/mean). If this value is under 0.5, coverage should be sufficiently even to allow confident base calling. If CV value is between 0.5 and 1, manual checking of base calling is recommended. If the CV value is greater than 1, some regions of the genome may not have sufficient coverage to allow nucleotide variance analysis.

24 Use the coverage values to correlate the qPCR results with sequencing success, such that future sequencing can be performed on samples with an appropriate number of genome copies to allow greatest chance of success.

25 Where mapping is not possible or advisable (when no referenced genome exists, or a large amount of novel accessory DNA is suspected), perform *de novo* assembly using SPAdes⁴⁶ or an equivalent program designed for use with data obtained from a sample which has undergone MDA. If required, align the resulting contigs against a reference sequence using abacas to determine core genomics regions, and to identify possible accessory or contaminant-derived data.

TIMING

Steps 1-15, IMS of bacteria from complex samples: minimum 2.5h

Steps 16-17, MDA: 2.5 h

Steps 18-25, DNA Quantification and Genome Sequencing: several days up to weeks, depending on method used

? TROUBLESHOOTING

Troubleshooting advice can be found in Table 1.

1

2 Anticipated Results

3 This protocol should yield complete bacterial genome sequences, derived through mapping
4 and variant nucleotide calling or *de novo* assembly. The ultimate success of the protocol can
5 only be assessed after analysis of the read data, although qPCR results can give an interim
6 indication as to the likely outcome. Sequence reads mapping to the target bacterium may
7 constitute only 9% of the output per sample, with the remaining reads representing
8 contaminating DNA, either human or microbial (Table 2). This is due to the fact that IMS is
9 used for enrichment of the target organism through depletion of contaminating DNA as
10 opposed to absolute purification of the target. The additional sequencing read data may
11 therefore provide information on the other organisms in the sample under analysis.

12 The success rate of the protocol is likely to depend on the nature and the load of target
13 organism in the input sample. Indeed repetition of the protocol on additional sample
14 material is unlikely to yield improved results and is not specifically recommended. Using *C.*
15 *trachomatis*, we achieved a success rate of 15-30% from urogenital samples and archived
16 diagnostic samples⁸. Rates were considerably higher for more recent, viable samples
17 archived after growth (100%)⁸. Using IMS-MDA, we have also been able to generate
18 complete genomes from older samples (up to 30 years old), which had been through several
19 rounds of freeze thawing and were no longer viable, with a 7% success rate (unpublished
20 data). For samples which do not produce complete genome sequences, partial sequence
21 data may still provide useful information for bioinformatic analysis.

22

23 Author contributions statements

24 N.R.T, P.M., I.N.C., H.S.-S. and S.H. conceived the project.

25 H.S.-S. designed and performed the experiments, researched and optimised the protocol and wrote
26 the manuscript.

27 N.R.T. supervised the development of the protocol, analysed the data and wrote the manuscript.

28 S.R.H. performed key data analysis, helped in the development of the protocol and wrote the
29 manuscript.

30 P.S. and S.P. contributed valuable information to the development of the protocol.

31 S.P. provided discarded clinical samples.

32 P.M., M.U., S.P., P.S., I.N.C. and J.P. provided insight into the experimental design and progress of
33 the protocol, and on the results and implications of the work.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35

Acknowledgments

This work was supported by the Wellcome Trust [grant number 098051]. The authors thank Louise Ellison for useful comments on the manuscript.

Competing financial interests

The authors declare that they have no competing financial interests.

Figure Legend

Figure 1. Schematic of the Immunomagnetic Separation and Multiple Displacement Amplification (IMS-MDA) procedure. In addition to the procedure shown, subsequent quantification, high throughput sequencing and data analysis are described in steps 18-25.

References

1. Mutreja, A. et al. Evidence for several waves of global transmission in the seventh cholera pandemic. *Nature* **477**, 462-5 (2011).
2. Bryant, J.M. et al. Whole-genome sequencing to identify transmission of *Mycobacterium abscessus* between patients with cystic fibrosis: a retrospective cohort study. *Lancet* **381**, 1551-60 (2013).
3. Harrison, E.M. et al. Whole genome sequencing identifies zoonotic transmission of MRSA isolates with the novel *mecA* homologue *mecC*. *EMBO Mol Med* **5**, 509-15 (2013).
4. Croucher, N.J., Harris, S.R., Grad, Y.H. & Hanage, W.P. Bacterial genomes in epidemiology - present and future. *Philos Trans R Soc Lond B Biol Sci* **368**, 20120202 (2013).
5. Holden, M.T. et al. A genomic portrait of the emergence, evolution, and global spread of a methicillin-resistant *Staphylococcus aureus* pandemic. *Genome Res* **23**, 653-64 (2013).
6. He, M. et al. Emergence and global spread of epidemic healthcare-associated *Clostridium difficile*. *Nat Genet* **45**, 109-13 (2013).
7. Eyre, D.W. et al. Detection of Mixed Infection from Bacterial Whole Genome Sequence Data Allows Assessment of Its Role in *Clostridium difficile* Transmission. *PLoS Comput Biol* **9**, e1003059 (2013).
8. Seth-Smith, H.M.B. et al. Whole-genome sequences of *Chlamydia trachomatis* directly from clinical samples without culture. *Genome Res* **23**, 855-866 (2013).
9. WHO, W.H.O. Prevalence and incidence of selected sexually transmitted infections. (World Health Organization, Geneva, 2011).

- 1 10. WHO, W.H.O. Global WHO alliance for the elimination of blinding trachoma by 2020. in
2 *Weekly epidemiological report* Vol. 17 (87) 161-8 (World Health Organisation, Geneva,
3 2012).
- 4 11. Dean, D. et al. Predicting phenotype and emerging strains among *Chlamydia trachomatis*
5 infections. *Emerg Infect Dis* **15**, 1385-94 (2009).
- 6 12. Klint, M. et al. High-resolution genotyping of *Chlamydia trachomatis* strains by multilocus
7 sequence analysis. *J Clin Microbiol* **45**, 1410-14 (2007).
- 8 13. Pannekoek, Y. et al. Multi locus sequence typing of Chlamydiales: clonal groupings within the
9 obligate intracellular bacteria *Chlamydia trachomatis*. *BMC Microbiol* **8**, 42 (2008).
- 10 14. Pedersen, L.N., Pødenphant, L. & Møller, J.K. Highly discriminative genotyping of *Chlamydia*
11 *trachomatis* using omp1 and a set of variable number tandem repeats. *Clin Microbiol Infect*
12 **14**, 644-652 (2008).
- 13 15. Jeffrey, B.M. et al. Genome sequencing of recent clinical *Chlamydia trachomatis* strains
14 identifies loci associated with tissue tropism and regions of apparent recombination. *Infect*
15 *Immun* **78**, 2544-53 (2010).
- 16 16. Harris, S.R. et al. Whole-genome analysis of diverse *Chlamydia trachomatis* strains identifies
17 phylogenetic relationships masked by current clinical typing. *Nat Genet* **44**, 364-6 (2012).
- 18 17. Joseph, S.J., Didelot, X., Gandhi, K., Dean, D. & Read, T.D. Interplay of recombination and
19 selection in the genomes of *Chlamydia trachomatis*. *Biol Direct* **6**, 28 (2011).
- 20 18. Andersson, P., Klein, M., Lilliebridge, R.A. & Giffard, P.M. Sequences of multiple bacterial
21 genomes and a *Chlamydia trachomatis* genotype from direct sequencing of DNA derived
22 from a vaginal swab diagnostic specimen. *Clin Microbiol Infect* **19**, E405-8 (2013).
- 23 19. Melnikov, A. et al. Hybrid selection for sequencing pathogen genomes from clinical samples.
24 *Genome Biol* **12**, R73 (2011).
- 25 20. Depledge, D.P. et al. Specific capture and whole-genome sequencing of viruses from clinical
26 samples. *PLoS One* **6**, e27805 (2011).
- 27 21. Bos, K.I. et al. A draft genome of *Yersinia pestis* from victims of the Black Death. *Nature* **478**,
28 506-10 (2011).
- 29 22. Kent, B.N. et al. Complete bacteriophage transfer in a bacterial endosymbiont (*Wolbachia*)
30 determined by targeted genome capture. *Genome Biol Evol* **3**, 209-18 (2011).
- 31 23. Geniez, S. et al. Targeted genome enrichment for efficient purification of endosymbiont DNA
32 from host DNA. *Symbiosis* **58**, 201-207 (2012).
- 33 24. Lasken, R.S. Genomic sequencing of uncultured microorganisms from single cells. *Nature Rev*
34 *Microbiol* **10**, 631-40 (2012).
- 35 25. Yilmaz, S. & Singh, A.K. Single cell genome sequencing. *Curr Op Biotechnol* **23**, 437-43 (2012).
- 36 26. Fitzsimons, M.S. et al. Nearly finished genomes produced using gel microdroplet culturing
37 reveal substantial intraspecies genomic diversity within the human microbiome. *Genome*
38 *Res* **23**, 878-88 (2013).
- 39 27. McLean, J.S. et al. Candidate phylum TM6 genome recovered from a hospital sink biofilm
40 provides genomic insights into this uncultivated phylum. *Proc Natl Acad Sci USA* **110**, E2390-
41 E2399 (2013).
- 42 28. Putman, T.E., Suchland, R.J., Ivanovitch, J.D. & Rockey, D.D. Culture-independent sequence
43 analysis of *Chlamydia trachomatis* in urogenital specimens identifies regions of
44 recombination and in-patient sequence mutations. *Microbiology* **10.1099/mic.0.070029-0**(2013).
- 45 29. Skjerve, E., Rørvik, L.M. & Olsvik, O. Detection of *Listeria monocytogenes* in foods by
46 immunomagnetic separation. *Appl Environ Microbiol* **56**, 3478-81 (1990).
- 47 30. Niesters, H.G.M., Quint, W.G.V. & Herbrink, P. Antigen capture directed polymerase chain
48 reaction for detection of *Chlamydia trachomatis*. in *Techniques in diagnostic pathology* Vol.
49 2 (eds. Bullock, G., Van Velzen, D. & Warhol, M.J.) (Academic Press Limited, 1991).
- 50

31. Fratamico, P.M., Schultz, F.J. & Buchanan, R.L. Rapid isolation of *Escherichia coli* O157:H7 from enrichment cultures of foods using an immunomagnetic separation method. *Food Microbiol* **9**, 105-13 (1992).
32. Hedrum, A., Lundeborg, J., Pålsson, C. & Uhlén, M. Immunomagnetic recovery of *Chlamydia trachomatis* from urine with subsequent colorimetric DNA detection. *Genome Res* **2**, 167-71 (1992).
33. Dean, F.B., Nelson, J.R., Giesler, T.L. & Lasken, R.S. Rapid amplification of plasmid and phage DNA using Phi29 DNA polymerase and multiply-primed rolling circle amplification. *Genome Res* **11**, 1095-9 (2001).
34. Hosono, S. et al. Unbiased whole-genome amplification directly from clinical samples. *Genome Res* **13**, 954-64 (2003).
35. Jalal, H. et al. Development and validation of a rotor-gene real-time PCR assay for detection, identification, and quantification of *Chlamydia trachomatis* in a single reaction. *J Clin Microbiol* **44**, 206-13 (2006).
36. Su, H., Watkins, N.G., Zhang, Y.X. & Caldwell, H.D. *Chlamydia trachomatis*-host cell interactions: role of the chlamydial major outer membrane protein as an adhesin. *Infect Immun* **58**, 1017-25 (1990).
37. Thornley, M.J. et al. Characterization of a monoclonal antibody to the group antigen of *Chlamydia* spp. and its use for antigen detection by reverse passive haemagglutination and indirect immunofluorescence. *FEMS Microbiol Lett* **17**, 45-9 (1983).
38. Thornley, M.J., Zamza, S.E., Byrne, M.D., Lusher, M. & Evans, R.T. Properties of monoclonal antibodies to the genus-specific antigen of *Chlamydia* and their use for antigen detection by reverse passive haemagglutination. *J Gen Microbiol* **131**, 7-15 (1985).
39. Zong, C., Lu, S., Chapman, A.R. & Xie, X.S. Genome-Wide Detection of Single-Nucleotide and Copy-Number Variations of a Single Human Cell. *Science* **338**, 1622-6 (2012).
40. Li, R., Li, Y., Kristiansen, K. & Wang, J. SOAP: short oligonucleotide alignment program. *Bioinformatics* **24**, 713-4 (2008).
41. Li, H. & Durbin, R. Fast and accurate short read alignment with Burrows-Wheeler transform. *Bioinformatics* **25**, 1754-60 (2009).
42. Carlson, J.H., Porcella, S.F., McClarty, G. & Caldwell, H.D. Comparative genomic analysis of *Chlamydia trachomatis* oculotropic and genitotropic strains. *Infect Immun* **73**, 6407-18 (2005).
43. Thomson, N.R. et al. *Chlamydia trachomatis*: genome sequence analysis of lymphogranuloma venereum isolates. *Genome Res* **18**, 161-71 (2008).
44. Stephens, R.S. et al. Genome sequence of an obligate intracellular pathogen of humans: *Chlamydia trachomatis*. *Science* **282**, 754-9 (1998).
45. Unemo, M. et al. The Swedish new variant of *Chlamydia trachomatis*: Genome sequence, morphology, cell tropism and phenotypic characterization. *Microbiology* **156**, 1394-1404 (2010).
46. Bankevich, A. et al. SPAdes: a new genome assembly algorithm and its applications to single-cell sequencing. *J Comput Biol* **19**, 455-77 (2012).
47. Chitsaz, H. et al. Efficient *de novo* assembly of single-cell bacterial genomes from short-read data sets. *Nature Biotechnol* **29**, 915-21 (2011).
48. Peng, Y., Leung, H.C.M., Yiu, S.M. & Chin, F.Y.L. IDBA-UD: A *de novo* assembler for single-cell and metagenomic sequencing data with highly uneven depth. *Bioinformatics* (2012).
49. Assefa, S., Keane, T.M., Otto, T.D., Newbold, C. & Berriman, M. ABACAS: Algorithm Based Automatic Contiguation of Assembled Sequences. *Bioinformatics* **25**, 1968-9 (2009).

1

2

1 **Table 1.** Troubleshooting table.

Step	Problem	Possible reason	Solution
2 & 13	Loss of beads in the reaction tube	Clumping of beads leading to poor retention on the magnet.	Ensure that Tween is added to PBST. Alternative buffers can also be used (see Dynabeads® manufacturer's instructions). Take care when removing buffers and ensure visually that the pellet of beads remains in the tube.
		Beads may be lost in the tube or pipette tip	Pulse spin tube if buffer is splashed up the side of tube. Avoid pipetting to mix beads: vortex or flick the tube.
10	Non-mixing of sample and bead mixture during incubation	The sample and bead mixture exist as separate droplets in the tube and do not come into contact during the incubation	Ensure that the bead-containing buffer is pipetted directly into the sample. If the sample is viscous, addition of further PBST may be advisable to a maximum of 500 µl.
21	Inefficient qPCR reaction	Primers / probe concentrations and cycling conditions inappropriate for assay	Use concentrations and conditions recommended for the specific assay. Use the standard dilutions to test a variety of conditions according to the manufacturer's instructions.
21	Insufficient amplification of target DNA	An inappropriate primary antibody may have been used, or insufficient quantities of it, in step 4.	This can be determined by retaining the supernatant in step 8 and performing qPCR to determine whether unbound target bacterium remains in the supernatant. Increased amounts of antibody or an alternative antibody can then be tested.
		Too much contaminating DNA or MDA reaction inhibitors remain in the samples	Perform two or more additional washes in PBST at step 14. A balance between removal of contamination and loss of or damage to the material of interest must be found.
		Excessive PBST buffer remains after IMS and is transferred to the MDA reaction	Ensure that as much PBST as possible is removed after the final wash (step 15). Ideally 1µl should remain, although MDA has been found to be effective even when 5µl buffer is added to the sample buffer with the beads.
		Poor sample quality	Select input (clinical) samples which are better preserved / of higher quality, or use a higher sample volume.

Zlotorynski, Eytan 30.8.13 10:23

Kommentar [4]: Au, OK or do you prefer WGA?

Helena 30.8.13 12:48

Kommentar [5]: Fine by me!