



**University of
Zurich**^{UZH}

**Zurich Open Repository and
Archive**

University of Zurich
University Library
Strickhofstrasse 39
CH-8057 Zurich
www.zora.uzh.ch

Year: 2016

Diagnostic accuracy of a smartphone electrocardiograph in dogs: Comparison with standard 6-lead electrocardiography

Vezzosi, T ; Buralli, C ; Marchesotti, F ; Porporato, F ; Tognetti, R ; Zini, E ; Domenech, O

Abstract: The diagnostic accuracy of a smartphone electrocardiograph (ECG) in evaluating heart rhythm and ECG measurements was evaluated in 166 dogs. A standard 6-lead ECG was acquired for 1 min in each dog. A smartphone ECG tracing was simultaneously recorded using a single-lead bipolar ECG recorder. All ECGs were reviewed by one blinded operator, who judged if tracings were acceptable for interpretation and assigned an electrocardiographic diagnosis. Agreement between smartphone and standard ECG in the interpretation of tracings was evaluated. Sensitivity and specificity for the detection of arrhythmia were calculated for the smartphone ECG. Smartphone ECG tracings were interpretable in 162/166 (97.6%) tracings. A perfect agreement between the smartphone and standard ECG was found in detecting bradycardia, tachycardia, ectopic beats and atrioventricular blocks. A very good agreement was found in detecting sinus rhythm versus non-sinus rhythm (100% sensitivity and 97.9% specificity). The smartphone ECG provided tracings that were adequate for analysis in most dogs, with an accurate assessment of heart rate, rhythm and common arrhythmias. The smartphone ECG represents an additional tool in the diagnosis of arrhythmias in dogs, but is not a substitute for a 6-lead ECG. Arrhythmias identified by the smartphone ECG should be followed up with a standard ECG before making clinical decisions.

DOI: <https://doi.org/10.1016/j.tvjl.2016.06.013>

Posted at the Zurich Open Repository and Archive, University of Zurich

ZORA URL: <https://doi.org/10.5167/uzh-128086>

Journal Article

Accepted Version



The following work is licensed under a Creative Commons: Attribution-NonCommercial-NoDerivatives 4.0 International (CC BY-NC-ND 4.0) License.

Originally published at:

Vezzosi, T; Buralli, C; Marchesotti, F; Porporato, F; Tognetti, R; Zini, E; Domenech, O (2016). Diagnostic accuracy of a smartphone electrocardiograph in dogs: Comparison with standard 6-lead electrocardiography. *Veterinary Journal*, 216:33-37.

DOI: <https://doi.org/10.1016/j.tvjl.2016.06.013>

**RELIABILITY OF A SMARTPHONE ELECTROCARDIOGRAM DEVICE IN
DOGS: COMPARISON WITH STANDARD 6-LEAD ELECTROCARDIOGRAPHY**

Tommaso Vezzosi, DVM; Carlotta Buralli, DVM; Federica Marchesotti, DVM; Federico Porporato, DVM; Rosalba Tognetti DVM, PhD; Eric Zini, DVM, PD, PhD; Oriol Domenech, DVM, MS.

From the Department of Veterinary Science (Vezzosi, Buralli, Tognetti), University of Pisa, Pisa, Italy; from the Department of Cardiology (Marchesotti, Porporato, Zini, Domenech), Istituto Veterinario di Novara, Novara, Italy; from the Clinic for Small Animal Internal Medicine, Vetsuisse Faculty, University of Zurich, Zurich, Switzerland (Zini); and from the Department of Animal Medicine, Production and Health, University of Padova, Legnaro, Italy (Zini).

The study was not supported by a grant.

Correspondence to: Dr. Vezzosi, tommaso.vezzosi@for.unipi.it

Abbreviations: HR, heart rate; App HR, heart rate automatically measured by the smartphone application; AVB, atrioventricular block; QTc, corrected QT interval.

ABSTRACT

Objective—To assess the reliability of a smartphone electrocardiogram (ECG) device in evaluating heart rhythm and ECG measurements in dogs.

Design—Prospective, multicenter, single-blind study.

Animals—166 client-owned dogs.

Procedures—A standard 6-lead ECG was acquired for 1 minute in each dog. A smartphone ECG tracing was simultaneously recorded using a single-lead bipolar ECG recorder. All ECGs were reviewed by one blinded operator, who judged if tracings were acceptable for interpretation and assigned an electrocardiographic diagnosis. Agreement between smartphone and standard ECG in the interpretation of tracings was evaluated. Sensitivity and specificity for the detection of arrhythmia were calculated with the smartphone ECG.

Results—Smartphone ECG tracings were interpretable in 162 out of 166 (97.6%) tracings. A perfect agreement between the smartphone and standard ECG was found in detecting bradycardia, tachycardia, ectopic beats and atrioventricular blocks. A very good agreement was found in detecting arrhythmias, with a 100% sensitivity and 97.9% specificity.

Conclusions and clinical relevance—The smartphone ECG provides tracings that are adequate for analysis in most dogs with a reliable assessment of heart rate, heart rhythm, atrioventricular blocks or ectopic beats. Remote smartphone ECG assessment can support the use of standard ECG in the management of dogs with arrhythmias.

Introduction

Many arrhythmias have paroxysmal presentation, while others require frequent monitoring due to the risk of progression. In these settings, serial electrocardiograms (ECG) are crucial for correct diagnosis and management. Clinical electrocardiography has thus undergone a continuous technological evolution since its invention by Willem Einthoven in 1903, leading to the development of Holter monitoring, telemetry systems and loop recorders.¹

One of the latest innovations is the 1-lead ECG recorded by smartphone devices using specific applications.²⁻⁴ Most smartphone ECG devices automatically digitize the electrocardiographic tracings. Some instruments enable PDF tracings to be uploaded to a cloud storage service or to send PDF files via email for remote interpretation. In human medicine there are many studies highlighting the good accuracy of smartphone ECG tracings to measure the heart rate (HR) and evaluate heart rhythm.⁵⁻⁷

Other studies have shown the good performances of smartphone ECG devices for QT interval assessment in patients receiving antiarrhythmic therapy,^{8,9} for diagnosing supraventricular tachycardia in pediatrics,¹⁰⁻¹² for detecting atrial fibrillation^{5,13-17} and for identifying signs associated with myocardial ischemia.^{18,19} Despite the considerable number of articles published in human medicine regarding the feasibility, accuracy and usefulness of smartphone ECG devices, to the best of our knowledge only one preliminary study has been performed comparing a smartphone ECG device to standardized ECG tracings in dogs.^a The aim of the present study was therefore to assess the use and reliability of a smartphone ECG to evaluate heart rhythm and ECG measurements in dogs.

Materials and methods

Animals

The study group included client-owned dogs that were referred to the Department of Veterinary Science of the University of Pisa or the Department of Cardiology of the Istituto Veterinario di Novara for a cardiologic consultation or assessment prior to anesthesia. The study was prospective, multicenter and single-blind. Dogs were recruited over a one-year period, from December 2014 to December 2015. Each case underwent a cardiac evaluation, including physical examination, standard 6-lead ECG and echocardiogram. The study protocol was reviewed and approved by the Institutional Welfare and Ethics Committee of the University of Pisa (permission number 39/2015). Written consent authorizing the participation of dogs in the study was obtained from each owner.

ECG acquisition and analysis

A standard 6-lead ECG^{b,c} was acquired for 1 minute in conscious, unsedated dogs that were positioned in right lateral recumbency. A smartphone ECG tracing was simultaneously recorded by three operators (TV, CB, FM) with an iPhone 4S^d using a single-lead bipolar ECG recorder^e and its application.^f The smartphone ECG was recorded placing it on the left precordial area of all the dogs. A cranio-caudal orientation of the smartphone case was used in each dog (Figure 1). In the short-haired dogs, a small amount of alcohol was placed on the left precordial area in order to obtain a good quality smartphone ECG signal. In the long-haired dogs, a small amount of alcohol was placed after shaving the left precordial area in order to acquire the same high-quality signal. The standard ECGs were stored digitally. Smartphone ECG recordings were automatically digitized by the device, sent via email and stored as a PDF. For each dog, ECG tracings obtained with the two methods were printed at a paper speed of 50 mm/s with a gain of 10 mm/mV. If QRS complexes had a very high or very low

amplitude, a 5 mm/mV or a 20 mm/mV gain was used, respectively. The last 30 s of each ECG tracing were analyzed. Dogs with a smartphone ECG trace lasting < 30 s were excluded from the study.

All ECG tracings were reviewed by a board-certified veterinary cardiologist (OD), in a blinded fashion, who judged whether the tracings were acceptable for interpretation. For all ECG tracings, the same operator evaluated the rhythm and performed ECG measurements. Measurements were achieved using the lead II of the standard ECG and using the only available lead of the smartphone ECG.

In each case, the following measurements were performed: mean HR (beats per minute, bpm); P wave amplitude (mv) and duration (ms); PQ interval duration (ms); R wave amplitude (mV); QRS complex duration (ms); QRS polarity; ST segment elevation or depression (mV); QT interval duration (ms); and corrected QT interval (QTc). The latter was calculated using the following formula: $QTc = \log 600 \times QT / \log RR$.²⁰ The mean HR was calculated by counting the number of QRS complexes in the last 30 s of the print out of each tracing. The result was multiplied by two in order to obtain the number of bpm. The QRS polarity of the smartphone ECG traces was compared with lead II of the standard ECG. Other measurements were achieved as previously described.²¹ Amplitude and duration measurements were calculated as the mean of three different beats. Finally, the mean HR calculated automatically by the smartphone application (App HR) was noted. Heart rate was classified as normal if between 70 and 160 bpm, bradycardia if < 70 bpm and tachycardia if > 160 bpm, as previously described.²¹

Statistical analysis

The analysis was performed only with paired ECG tracings that were considered acceptable for interpretation, as defined by the operator, and the standard ECG was set as the reference

method. Cohen's kappa (κ) test was used to calculate the agreement between the smartphone ECG and standard ECG for HR classification (normal, bradycardia or tachycardia), heart rhythm (sinus rhythm, atrial fibrillation, ventricular rhythm, supraventricular rhythm), atrioventricular blocks (AVB) (absent, first-degree AVB, second-degree AVB, third-degree AVB), premature complexes (absent, ventricular, supraventricular), ST segment (normal, elevated or depressed), polarity of QRS complex (positive, negative) and QTc interval (normal, long or short). The kappa coefficient was interpreted as follows: values ≤ 0.20 as no agreement, 0.21–0.40 as fair, 0.41–0.60 as moderate, 0.61–0.80 as good, 0.81–0.99 as very good, and 1 as perfect agreement.

The sensitivity, specificity, positive and negative predictive values of the smartphone ECG to detect arrhythmia were calculated. In addition, the median and range of differences between the standard ECG and smartphone ECG were calculated for the amplitude of the P and R waves, for the duration of the P wave, PQ interval, QRS complex and QTc interval. The Pearson or Spearman correlation coefficients were used to study correlations between HRs measured with the standard ECG and the smartphone ECG using values calculated either by the operator or automatically. The Shapiro-Wilk test was used to determine the normality of the datasets. Statistical analysis was performed with a commercial software.^g $P < 0.05$ was considered as significant.

Results

Animals and feasibility

A total of 166 dogs were enrolled in the study, of which 84 were males and 82 were females. The median age was 9 years (ranging between 0.3 to 17 years) and median body weight was 25 kg (55.1 lb) with a range of 2.1 to 75 kg (4.6 to 165.3 lb). The majority of dogs (71 out of 166, 43%) had cardiac diseases, both congenital or acquired; 32 dogs (19%) had neoplastic

diseases, 31 (18%) were in the intensive care unit because of renal, respiratory, gastro-intestinal or neurologic diseases, and 33 (20%) were healthy dogs evaluated for pre-anesthesia assessment prior to elective surgeries.

The blinded operator judged 162 of 166 (97.6%) of the smartphone ECG tracings as acceptable for interpretation (Figures 2, 3 and 4). In 4 cases (2.4%), the tracings were deemed as non-interpretable; all of tracings were recorded in dogs with a body weight < 10 kg.

Heart rate

According to the standard 6-lead ECG, 133 out of 162 dogs (82%) had a normal HR, 20 (12%) had tachycardia, and 9 (6%) had bradycardia. A perfect agreement ($\kappa=1$) between the smartphone and standard ECG was found in the classification of HR when it was manually measured on digitized tracings (Table 1). A strong positive correlation was found between the HR values obtained by both methods ($r^2 = 0.99$; $p < 0.0001$; Figure 5). Median paired differences between the HR manually measured on standard ECG and smartphone ECG was 0 bpm (-10, +20 bpm; Table 2 and Figure 6). A strong positive correlation was also found between the App HR values and those manually measured on standard ECG tracings ($r^2 = 0.923$; $p < 0.0001$; Figure 7). However, the App HR was less accurate than the manually measured HR on digitized smartphone ECG tracings ($\kappa=0.91$). In fact, in 103 out of 162 (63.6%) cases, the App HR underestimated the actual HR, with a median difference of -3 bpm; (-31, +20 bpm; Figure 8). However in only 4/162 (2.5%) cases, was there a misclassification of HR with the smartphone application. According to App HR, two dogs with tachycardia were classified as normal HR, one dog with normal HR was classified as bradycardia, and one dog with bradycardia was classified as normal HR. The greatest disagreement was found in a dog with severe bradycardia (40 bpm) secondary to a third-

degree AVB in which the App HR read the P waves as QRS complexes, thus erroneously yielding an HR of 140 bpm.

Heart rhythm

One hundred and forty-one dogs (87%) had sinus rhythm or sinus arrhythmia; 14 dogs (9%) had supraventricular arrhythmias; 7 dogs (4%) had ventricular rhythm or ventricular arrhythmias. Six dogs (4%) had different types of AVBs.

Very good agreement ($\kappa=0.94$) was found in the evaluation of the heart rhythm. Disagreement was found in only 3 out of 162 (1.9%) cases, in which the sinus rhythm was erroneously classified as atrial rhythm due to the negative polarity of the P waves (1 case) or as a slow atrial fibrillation due to non observable P waves (2 cases) on the smartphone ECG trace (Table 3). In 128 out of 149 (85.9%) cases, the smartphone ECG underestimated the amplitude of the P wave, with a median difference of -0.1 mV (-0.4; +0.1 mV). The analysis of the P wave duration showed a median difference between the two methods of 0 ms (-20, +0 ms).

QRS complex analysis

A good agreement ($\kappa=0.65$) was found in the polarity of the QRS complexes between the smartphone ECG and lead II of the standard 6-lead ECG (Figures 2, 3, 4). The same QRS polarity was found in 158 out of 162 (97.5%) cases. In 3 cases with positive polarity of the QRS complex in lead II, the smartphone tracing showed a negative QRS. In 1 case with negative polarity of the QRS complex in lead II, the smartphone tracing showed a positive QRS. The evaluation of the QRS duration showed a median difference of 0 ms (-20, +10 mV). Lastly, the smartphone ECG underestimated the amplitude of R wave in 121 out of 162 (74.7%), with a median difference of -0.5 mV (-2.1; +1), compared to the standard ECG.

188

189 *Ectopic beats*

190 A perfect agreement ($\kappa=1$) between the smartphone ECG and standard ECG was found in the
191 identification and classification of ectopic beats, including 16 cases with ventricular
192 premature complexes, 3 cases with supraventricular premature complexes and 4 cases with
193 both supraventricular and ventricular ectopic beats. In addition a perfect agreement was found
194 regarding the polarity of ventricular premature complexes on the smartphone ECG tracings
195 compared with lead II of the standard 6-lead ECG.

196

197 *Atrioventricular blocks*

198 A perfect agreement ($\kappa=1$) between the smartphone ECG and standard ECG was found in the
199 AVB diagnosis, including 2 cases with first-degree AVB, 1 with second-degree AVB and 3
200 cases with third-degree AVB. The PQ interval analysis using smartphone tracings was
201 reliable in comparison to the standard ECG, with a median difference of 0 ms (range -20, +20
202 ms).

203

204 *QT interval and ST segment*

205 In the interpretation of the ST segment and QT interval, the agreement between the
206 smartphone ECG and standard ECG was good ($\kappa=0.70$ and $\kappa=0.72$, respectively). In 2 dogs
207 with a normal ST segment on the standard ECG, the smartphone ECG erroneously indicated
208 ST depression in 1 case and ST elevation in the other. In 2 dogs the smartphone ECG did not
209 identify the ST depression, which was evident on the standard ECG. In our study, the
210 smartphone ECG underestimated the QT interval, with a median difference of -10 ms (range -
211 34, +19 ms) compared to the standard ECG. Five cases of long QT intervals were erroneously
212 classified as normal QT intervals using tracings recorded with the smartphone ECG.

However, in 7 cases with moderately to severely increased QT intervals ($QT_c > 260$ ms), the long QT was correctly identified.

Performance with arrhythmias

Considering all the arrhythmias taken together, the smartphone ECG had 100% sensitivity and 97.9% specificity in differentiating between sinus rhythm and arrhythmias, with a positive predictive value of 87.5%, and a negative predictive value of 100%.

Discussion

In our investigation the smartphone ECG was easily performed in all dogs and 96.7% of tracings were deemed as interpretable. These results are in line with findings in human medicine where smartphone ECG tracings were interpretable in 87-99% of patients.^{12,13,22} The few tracings judged as non-interpretable were all recorded in small breed dogs, where motion artifacts are common, which likely accounted for the fact that the tracings were not readable. The smartphone ECG was excellent in the HR evaluation in dogs, since a strong positive correlation was found between the HR values obtained by the smartphone device and the standard ECG. This is in accordance with a preliminary study in dogs, where instantaneous and average heart rates were identical in all cases where exact matches were possible in a comparison between smartphone ECG and reference ECG.^a

In our investigation, the greatest reliability was found when the HR was manually measured on digitized tracings. Conversely the App HR was less reliable, since lower agreement was found between the HR values obtained by the smartphone device and the standard ECG. As the QRS complexes on smartphone ECG tracings had a low amplitude in most dogs, we hypothesize that the App HR may underestimate the HR due to the fact that some QRS complexes are not correctly interpreted by the instrument. In a few dogs, the App HR was

totally unreliable. However, in only 1 case the disagreement was of a real clinical value: in a dog with severe bradycardia secondary to third-degree AVB, the App HR read the P waves as QRS complexes, thus erroneously resulting in a normal HR.

The smartphone ECG was very reliable in evaluating heart rhythm in dogs, as it showed 100% sensitivity and 98% specificity in the detection of arrhythmias. The recorded arrhythmias were also accurately identified using the smartphone ECG, in most instances. All cases of atrial fibrillation were correctly diagnosed, without false negatives. This result is similar to findings in humans where the sensitivity and specificity of the smartphone ECG in detecting atrial fibrillation were 94-100% and 90-97%, respectively.^{5,7,22} In humans, most false diagnoses of atrial fibrillation are due to small voltage P waves. Our results showed that the smartphone ECG underestimates the amplitude of the P wave. Despite this, the P waves remained clearly visible in the majority of dogs with sinus rhythm. In a few cases, however, the P waves were difficult to recognize and occasionally it was hard to differentiate between sinus rhythm and atrial fibrillation. Consequently, 2 out of 141 cases of sinus rhythm were incorrectly classified as atrial fibrillation. In a small breed dog, the P waves had negative polarity on the smartphone ECG leading to the incorrect diagnosis of an atrial ectopic rhythm.

A preliminary study in cats recommended positioning the smartphone case parallel to the long axis of the heart, with a more base-apex orientation in comparison to our cranio-caudal orientation.^h It might be that in some small breed dogs, the orientation of the smartphone case should be individually adjusted to correctly visualize the P waves. Atrial fibrillation is common in dogs with severe cardiac disease and increases the risk of cardiac-related death in those with myxomatous mitral valve degeneration and dilated cardiomyopathy.^{23,24} Likewise, in humans, atrial fibrillation increases the chance of morbidity or mortality, and recent studies have highlighted the utility of the smartphone ECG in screening for this arrhythmia.^{5,7,13-}

262 ^{16,25,26} A significant number of occult atrial fibrillations have been detected by physicians,
263 pharmacists or directly by patients using an ECG recorded with a smartphone.

264 Early diagnosis of atrial fibrillation is difficult in dogs. Our results show that the smartphone
265 ECG may become a promising tool for frequent monitoring of dogs predisposed to atrial
266 fibrillation. It could also be beneficial for dogs with atrial fibrillation that receive drugs to
267 control HR. Holter monitoring is an essential tool for evaluating HR and in treating atrial
268 fibrillation in dogs. However, 24-hour Holter monitoring is expensive and necessitates the
269 owner's compliance, hence its use may not always be practical. In the light of its ease and
270 cost effectiveness, the smartphone ECG could represent a complementary tool for HR
271 evaluation at home in dogs with atrial fibrillation.

272 The smartphone ECG showed a good reliability in the analysis of the QRS complex, in
273 assessing both duration and polarity. In most dogs, QRS complexes displayed the same
274 polarity on smartphone tracings and lead II of the 6-lead ECG, with a similar polarity in all
275 cases of ventricular ectopic beats. In comparison to the standard ECG however, the
276 smartphone device underestimated the R wave amplitude. In many dogs it was necessary to
277 use the highest calibration setting (20 mm/mV) to optimally visualize the electrocardiographic
278 waves. Therefore, smartphone tracings should not be used to assess the amplitude of ECG
279 waves as a substitute for standard electrocardiograms.

280 The smartphone ECG was highly reliable in the identification of ectopic beats. Ventricular
281 premature complexes, accelerated idioventricular rhythms and ventricular tachycardias were
282 easily identified in all cases with the smartphone ECG. One recent investigation used it as the
283 sole electrocardiographic method in the identification of ventricular premature complexes in
284 the screening of Doberman Pinschers for occult dilated cardiomyopathy.¹ It could thus be
285 useful in screening or monitoring dogs with cardiomyopathies associated with ventricular
286 arrhythmias.

With regard to the reliability of the smartphone ECG for AVBs, a good agreement with the standard ECG was found both in the evaluation of the PQ interval and in the identification of the type of block. One study in humans described a higher percentage of false positives and negatives during the evaluation of AVBs compared to our results.⁷ The authors reported motion artifacts (arm movement, muscle tension and tremor) as the main difficulties in AVB evaluation. None of the smartphone ECG tracings recorded motion artifacts that led to misdiagnosing AVBs. Thus, the agreement between devices was perfect, suggesting that the smartphone ECG can be helpful in the interpretation of AVBs in dogs.

Lastly, the reliability of the smartphone device in the assessment of the QT interval and ST segment was evaluated. The smartphone ECG underestimated the QT interval in most instances, leading to false negative results in the diagnosis of the long QT. However, in all dogs with an increased QT ($QT_c > 260$ ms), there was a complete agreement. In one study performed on humans, the smartphone ECG was highly reliable in the analysis of the QT interval.⁸ In our group, the morphology of the T wave based on the smartphone ECG was not always similar to that reported with the standard ECG. This finding, along with the frequent biphasic morphology of the T wave, makes the evaluation of the QT interval less accurate in dogs. Misclassifications were also found in the ST segment interpretation. Taken together, these findings suggest that the smartphone ECG may only be partially reliable for QT interval and ST segment evaluation in dogs.

Our investigation has some limitations. First, the study group was large but the number of dogs with arrhythmias was relatively low. A larger number of rhythm disturbances might have revealed a lower reliability of the smartphone ECG. However, most common types of canine arrhythmias were included in our study and in all these cases the smartphone ECG was consistent in diagnosing the arrhythmia. Second, during the recording of the smartphone ECG, all dogs were positioned in right lateral recumbency and it is not known if there would

have been similar results if the left side or standing positioning was used instead. Furthermore, we recorded 1 minute with both ECG methods, but only the last 30 s were used for the manual HR assessment. As the App HR provided an estimate of the mean HR based on a 1-minute registration, we cannot rule out that the interpretation of the results was partly biased.

In conclusion, the smartphone ECG can rapidly and simply record a single-lead ECG of good diagnostic quality in dogs. Tracing analysis performed by cardiologists reliably evaluated HR, heart rhythm, AVBs and ectopic beats. The smartphone ECG does not substitute the 6-lead ECG or Holter monitoring but does represent an additional tool in the management of dogs with arrhythmias or in monitoring dogs at risk for heart rhythm disturbances. Further studies are needed to assess the diagnostic value of the smartphone ECG recorded by owners in a home setting.

Footnotes

^a Kraus MS, Rishniw MR, Gelzer AR, et al. Comparison of the AliveCor® ECG device for the iPhone with a reference standard electrocardiogram. ACVIM forum 2013. Seattle (WA), USA.

^b Elan 1100 ECG system, Cardioline, et medical devices SpA, Milano, Italy.

^c MAC 800 ECG system, GE Healthcare, Milano, Italy.

^d AliveCor Veterinary Heart Monitor, AliveCor, Inc., San Francisco, CA, USA.

^e AliveECG Vet, AliveCor, Inc., San Francisco, CA, USA.

^f Apple Inc, Cupertino, CA, USA.

^g GraphPad Prism 5, San Diego, CA, USA.

^h Stromberg V, Kvart C. iPhone ECG and Hypertrophic Cardiomyopathy (HCM) in Cats. Epsilon Archive for Student Projects. Library of the Swedish University of Agricultural Sciences. Uppsala. 2015. Online link: <http://stud.epsilon.slu.se/7583> (last access on 06-02-2016).

ⁱ Gordon S, Estrada AH, Braz-Ruivo L, et al.. Evaluation of NTproBNP, High Sensitivity Troponin I and PDK4 for the Detection of Occult DCM: A Prospective Study in 449 Doberman Pinschers. 25th ECVIM-CA Congress 2015. Lisbon. Portugal.

Conflict of Interest

The authors disclose no conflict of interest.

References

1. Kennedy HL. The evolution of ambulatory ECG monitoring. *Prog Cardiovasc Dis* 2013; 56:127-132.
2. Bruining N, Caiani E, Chronaki C, et al. Task Force of the e-Cardiology Working. Acquisition and analysis of cardiovascular signals on smartphones: potential, pitfalls and perspectives: by the Task Force of the e-Cardiology Working Group of European Society of Cardiology. *Eur J Prev Cardiol* 2014;21(2 Suppl):4-13.
3. Walsh JA 3rd, Topol EJ, Steinhubl SR. Novel wireless devices for cardiac monitoring. *Circulation* 2014;130:573-581.
4. Baquero GA, Banchs JE, Ahmed S, et al. Surface 12 lead electrocardiogram recordings using smart phone technology. *J Electrocardiol* 2015;48:1-7.
5. Lau JK, Lowres N, Neubeck L, et al. iPhone ECG application for community screening to detect silent atrial fibrillation: a novel technology to prevent stroke. *Int J Cardiol* 2013;165:193-194.

6. Ho CL, Fu YC, Lin MC, et al. Smartphone applications (apps) for heart rate measurement in children: comparison with electrocardiography monitor. *Pediatr Cardiol* 2014;35:726-731.
7. Haberman ZC, Jahn RT, Bose R, et al. Wireless Smartphone ECG Enables Large-Scale Screening in Diverse Populations. *J Cardiovasc Electrophysiol* 2015;26:520-526.
8. Chung EH, Guise KD. QTC intervals can be assessed with the AliveCor heart monitor in patients on dofetilide for atrial fibrillation. *J Electrocardiol* 2015;48:8-9.
9. Muhlestein JB. QTC intervals can be assessed with the AliveCor heart monitor in patients on dofetilide for atrial fibrillation. *J Electrocardiol* 2015;48:10-11.
10. Wackel P, Beerman L, West L, et al. Tachycardia detection using smartphone applications in pediatric patients. *J Pediatr* 2014;164:1133-1135.
11. Ferdman DJ, Liberman L, Silver ES. A Smartphone Application to Diagnose the Mechanism of Pediatric Supraventricular Tachycardia. *Pediatr Cardiol* 2015;36:1452-1457.
12. Nguyen HH, Van Hare GF, Rudokas M, et al. SPEAR Trial: Smartphone Pediatric ElectroCARDiogram Trial. *PLoS One* 2015;10:e0136256.
13. Saxon LA. Ubiquitous wireless ECG recording: a powerful tool physicians should embrace. *J Cardiovasc Electrophysiol* 2013;24:480-483.
14. Lee J, Reyes BA, McManus DD, et al. Atrial fibrillation detection using an iPhone 4S. *IEEE Trans Biomed Eng* 2013;60:203-206.
15. McManus DD, Lee J, Maitas O, et al. A novel application for the detection of an irregular pulse using an iPhone 4S in patients with atrial fibrillation. *Heart Rhythm* 2013;10:315-319.
16. Orchard J, Freedman SB, Lowres N, et al. iPhone ECG screening by practice nurses

and receptionists for atrial fibrillation in general practice: the GP-SEARCH qualitative pilot study. *Aust Fam Physician* 2014; 43:315-319.

17. Lowres N, Krass I, Neubeck L, et al. Atrial fibrillation screening in pharmacies using an iPhone ECG: a qualitative review of implementation. *Int J Clin Pharm* 2015;37:1111-1120.

18. Wong CK. iPhone ECG monitoring: the gateway to the new paradigm of STEMI therapy. *Int J Cardiol* 2013;168:2897-2898.

19. Muhlestein JB, Le V, Albert D, et al. Smartphone ECG for evaluation of STEMI: results of the ST LEUIS Pilot Study. *J Electrocardiol* 2015;48:249-259.

20. Matsunaga T, Mitsui T, Harada T, et al. QT corrected for heart rate and relation between QT and RR intervals in beagle dogs. *J Pharmacol Toxicol Methods* 1997;38:201-209.

21. Kittleson MD, Kienle RD. Small animal cardiovascular medicine. Mosby Inc, 1998; 85.

22. Tarakji KG, Wazni OM, Callahan T, et al. Using a novel wireless system for monitoring patients after the atrial fibrillation ablation procedure: the iTransmit study. *Heart Rhythm*. 2015;12:554-559.

23. Calvert CA, Pickus CW, Jacobs GJ, et al. Signalment, survival, and prognostic factors in Doberman pinschers with end-stage cardiomyopathy. *J Vet Intern Med* 1997;116:323-326.

24. Jung SW, Sun W, Griffiths LG, et al. Atrial Fibrillation as a Prognostic Indicator in Medium to Large-Sized Dogs with Myxomatous Mitral Valvular Degeneration and Congestive Heart Failure. *J Vet Intern Med* 2016;30:51-57.

25. Peritz DC, Howard A, Ciocca M, et al. Smartphone ECG aids real time diagnosis of palpitations in the competitive college athlete. *J Electrocardiol* 2015;48:896-899.

411 26. Lowres N, Freedman SB, Gallagher R, et al. Identifying postoperative atrial
412 fibrillation in cardiac surgical patients posthospital discharge, using iPhone ECG: a
413 study protocol. *BMJ Open* 2015; 5:e006849.

414

Tables

Table 1. Agreement (κ) between smartphone ECG and standard 6-lead ECG.

Type of analysis	κ	(95% CI)	Agreement
Manual HR	1		Perfect
App HR	0.91	(0.81-0.99)	Very good
Heart rhythm	0.94	(0.86-1)	Very good
AVBs	1		Perfect
Ectopic beats	1		Perfect
ST interval	0.70	(0.43-0.98)	Good
QTc interval	0.72	(0.49-0.95)	Good
QRS polarity	0.65	(0.34 to 0.97)	Good

CI = confidence interval; Manual HR = HR manually measured on printed ECG tracings; App HR = HR automatically measured by smartphone application; AVBs = atrioventricular blocks; QTc = corrected QT interval.

Table 2. Differences between smartphone ECG and standard ECG in the evaluation of electrocardiographic parameters.

Parameter	Difference	Range
Manual HR (bpm)	0	-10; +20
App HR (bpm)	-3	-31; +20
P (ms)	0	-20; +0
P (mV)	-0,1	-0,4; +0,1
PQ (ms)	0	-20; +20
QRS (ms)	0	-20; +10
R (mV)	-0,5	-2,1; +1
QTc (ms)	-10	-34; +19

Manual HR = HR manually measured on printed ECG tracings; App HR = HR automatically measured by smartphone application. Median difference and range are reported.

427 **Table 3.** Agreement between smartphone ECG and standard 6-lead ECG in heart rhythm
 428 identification in 162 dogs.

Standard ECG	Smartphone ECG				Total
	S	AF	SV	V	
S	138	2	1	0	141
AF	0	12	0	0	12
SV	0	0	2	0	2
V	0	0	0	7	7
Total	138	14	3	7	162

429 S, sinus rhythm; AF, atrial fibrillation; SV, supraventricular rhythm; V, ventricular rhythm.

430

Figure legends

Figure 1. Cranio-caudal orientation of the smartphone in a dog. The camera side of the smartphone was located caudally.

Figure 2. Sinus rhythm with standard ECG (A) and with smartphone ECG (B) in the same dog. Paper speed = 50 mm/s; 10 mm = 1 mV.

Figure 3. Atrial fibrillation with standard ECG (A) and with smartphone ECG (B) in the same dog. Paper speed = 50 mm/s; 10 mm = 1 mV.

Figure 4. Third-degree AVB with standard ECG (A) and with smartphone ECG (B) in the same dog. Paper speed = 50 mm/s; 5 mm = 1 mV.

Figure 5. Pearson test showing a strong positive correlation between the HR values manually measured on standard ECGs and smartphone ECG tracings ($r^2 = 0.99$; $p < 0.0001$).

Figure 6. Bland-Altman plot showing differences between HR values manually measured on standard ECG and smartphone ECG tracings.

Figure 7. Pearson test showing a strong positive correlation between the HR values manually measured on standard ECGs and HR values produced by the smartphone application ($r^2 = 0.92$; $p < 0.0001$).

455 **Figure 8.** Bland-Altman plot showing differences between the HR values manually measured
456 on standard ECGs and HR values produced by the smartphone application.