



**University of  
Zurich<sup>UZH</sup>**

**Zurich Open Repository and  
Archive**

University of Zurich  
University Library  
Strickhofstrasse 39  
CH-8057 Zurich  
[www.zora.uzh.ch](http://www.zora.uzh.ch)

---

Year: 2016

---

## **Performance of the Assurance GDS® assay for the detection of *L. monocytogenes* in pure cultures and spiked food samples**

Althaus, Denise ; Zweifel, Claudio ; Johler, Sophia ; Stephan, Roger

**Abstract:** *Listeria monocytogenes* is a foodborne pathogen with significant impacts on public health and economy worldwide. Reliable and fast detection of *L. monocytogenes* is of major importance for both diagnostic laboratories and the food industry. The current study evaluated the performance of the Assurance GDS® assay for the detection of *L. monocytogenes* in pure cultures and spiked food samples. In the pure culture experiments, the Assurance GDS® assay for *Listeria monocytogenes* accurately detected the target strains of different serotypes and was correctly negative for a variety of other *Listeria* species. For reliable detection of *L. monocytogenes* in pure culture experiments, colony counts >10<sup>5</sup> cfu/ml were required, which emphasizes the need for an adequate enrichment step. The challenge test experiments (steak tartare, bologna type sausage, Gorgonzola cheese) using a one-broth enrichment strategy showed that the Assurance GDS® assay reliably detected *L. monocytogenes* after 16 h of enrichment in Half-Fraser broth, provided that spiking levels of the different matrices were 10<sup>2</sup> cfu/g. Depending of the food matrix, longer incubation times of 24 h or 48 h were required when the initial spiking level was <10<sup>2</sup> cfu/g, as to be expected in a proportion of naturally contaminated food products. Thus, the Assurance GDS® *Listeria monocytogenes* assay has proven to be a reliable and easy to handle, rapid test system for the specific detection of *L. monocytogenes*. This system is a suitable tool for generating microbiological results used for a “positive batch release”, especially for RTE foods with short shelf lives. However, longer enrichment times (24 h or 48 h) are required in a one-broth enrichment strategy, when the contamination level of the food matrix is low (<10<sup>2</sup> cfu/g).

DOI: <https://doi.org/10.4172/2476-2059.1000103>

Posted at the Zurich Open Repository and Archive, University of Zurich

ZORA URL: <https://doi.org/10.5167/uzh-129859>

Journal Article

Published Version



The following work is licensed under a Creative Commons: Attribution 4.0 International (CC BY 4.0) License.

Originally published at:

Althaus, Denise; Zweifel, Claudio; Johler, Sophia; Stephan, Roger (2016). Performance of the Assurance GDS® assay for the detection of *L. monocytogenes* in pure cultures and spiked food samples. *Journal of Food: Microbiology, Safety and Hygiene*, 1:103.

DOI: <https://doi.org/10.4172/2476-2059.1000103>



## Performance of the Assurance GDS® Assay for the Detection of *L. monocytogenes* in Pure Cultures and Spiked Food Samples

Denise Althaus, Claudio Zweifel, Sophia Johler and Roger Stephan\*

Institute for Food Safety and Hygiene, Vetsuisse Faculty University of Zurich, Zurich, Switzerland

\*Corresponding author: Roger Stephan, Institute for Food Safety and Hygiene, Vetsuisse Faculty University of Zurich, Winterthurerstrasse 272, 8057 Zurich, Switzerland, Tel: +41-44-635-8651; Fax: +41-44-635-8908; E-mail: stephanr@fsafety.uzh.ch

Received date: January 04, 2016; Accepted date: February 29, 2016; Published date: March 04, 2016

Copyright: © 2016 Althaus D, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

### Abstract

*Listeria monocytogenes* is a foodborne pathogen with significant impacts on public health and economy worldwide. Reliable and fast detection of *L. monocytogenes* is of major importance for both diagnostic laboratories and the food industry. The current study evaluated the performance of the Assurance GDS® assay for the detection of *L. monocytogenes* in pure cultures and spiked food samples. In the pure culture experiments, the Assurance GDS® assay for *Listeria monocytogenes* accurately detected the target strains of different serotypes and was correctly negative for a variety of other *Listeria* species. For reliable detection of *L. monocytogenes* in pure culture experiments, colony counts  $>10^5$  cfu/ml were required, which emphasizes the need for an adequate enrichment step. The challenge test experiments (steak tartare, bologna type sausage, Gorgonzola cheese) using a one-broth enrichment strategy showed that the Assurance GDS® assay reliably detected *L. monocytogenes* after 16 h of enrichment in Half-Fraser broth, provided that spiking levels of the different matrices were  $\geq 10^2$  cfu/g. Depending of the food matrix, longer incubation times of 24 h or 48 h were required when the initial spiking level was  $<10^2$  cfu/g, as to be expected in a proportion of naturally contaminated food products. Thus, the Assurance GDS® *Listeria monocytogenes* assay has proven to be a reliable and easy to handle, rapid test system for the specific detection of *L. monocytogenes*. This system is a suitable tool for generating microbiological results used for a "positive batch release", especially for RTE foods with short shelf lives. However, longer enrichment times (24 h or 48 h) are required in a one-broth enrichment strategy, when the contamination level of the food matrix is low ( $<10^2$  cfu/g).

**Keywords:** Assurance GDS® detection system; *Listeria monocytogenes*; Performance; Pure cultures; Spiked food samples; Enrichment times

### Introduction

*Listeria monocytogenes* is an important foodborne pathogen that has significant impacts on public health and economy worldwide. *L. monocytogenes* belongs to the genus *Listeria*, which actually includes 18 further species: *L. aquatica*, *L. booriae*, *L. cornellensis*, *L. denitrificans*, *L. fleischmannii*, *L. floridensis*, *L. grandensis*, *L. grayi*, *L. innocua*, *L. ivanovii*, *L. marthii*, *L. murrayi*, *L. newyorkensis*, *L. riparia*, *L. rocourtiae*, *L. seeligeri*, *L. weihenstephanensis*, and *L. welshimeri* (<http://www.bacterio.net>). *L. monocytogenes* has the potential to cause serious and life-threatening conditions (including septicemia, meningitis, meningoencephalitis, and abortion) in persons with reduced immunity [1]. In the European Union, a total of 2.161 confirmed human cases of listeriosis (notification rate of 0.52 cases per 100.000 population) were reported in 2014 [2]. Ready-to-eat (RTE) foods seem to cause the majority of human *L. monocytogenes* infections and RTE products have also been implicated in large-scale outbreaks [3-6]. *Listeria* spp. are widely distributed in the environment and certain strains may become established and persist in the processing environment [7-9].

Detection of *L. monocytogenes* traditionally involves culture methods including selective enrichment and plating, followed by the biochemical identification of presumptive *L. monocytogenes* colonies. Using culture-based techniques, as e.g. ISO/EN 11290-1 [10], it takes up to one week until the identification is completed. However, there is

a growing need for rapid tests generating results comparable to standard methods. Such rapid tests are of special importance for products with short shelf lives and (real-time) batch releases by food processing companies. Hence, several immunological and molecular biological methods have been developed [11]. Commercially available rapid molecular detection systems for *L. monocytogenes* or *Listeria* spp. include e.g. the Assurance GDS® assays, the GeneQuence® assay, the iQ-Check® kit, the Qualicon BAX® system, or the TaqMan® detection kit. The Assurance GDS® assay thereby combines the PCR approach with a preceding immunomagnetic separation (IMS) step [12,13]. The aim of the present study was (i) to determine the diagnostic specificity and sensitivity of the Assurance GDS® assay for *L. monocytogenes* and (ii) to evaluate the performance of this system for detection of *L. monocytogenes* in selected ready-to-eat food products using a one-broth enrichment strategy.

### Materials and Methods

#### Specificity of the assurance GDS® *L. monocytogenes* assay

To determine the specificity of the Assurance GDS® *L. monocytogenes* kit (Bio Control Systems, Bellevue, WA, USA), pure culture experiments were performed using a collection of seven *L. monocytogenes* strains (serotypes 1/2a, 1/2b, 1/2c, 3a, 3c, 4b, O-group 4 non-motile) and 13 strains of various other *Listeria* species (Table 1). Strains originated from the collection of the Institute for Food Safety and Hygiene (University of Zurich, Switzerland) or were obtained from the Leibniz Institute DSMZ (German Collection of Microorganisms and Cell Cultures, Braunschweig, Germany). Moreover, the serotype 3a

*L. monocytogenes* strain was obtained from the BC Centre for Disease Control (BCCDC, Vancouver, Canada). After growth of the 20 strains on sheep blood agar (overnight at 37°C; Oxoid, Pratteln, Switzerland), single colonies were inoculated into 10 ml of brain heart infusion broth (BHI; Oxoid) and incubated overnight (16 h) at 37°C. Colony counting (plate count agar, 24 h at 37°C; Oxoid) confirmed that this procedure yielded stationary phase cultures containing about 10<sup>9</sup> cfu/ml. Subsets of the incubated BHI broth cultures were tested using the Assurance GDS® *L. monocytogenes* (Bio Control Systems) according to the manufacturers' instructions.

Species	Sero-type	Strain designation	Assurance GDS® <i>L. monocytogenes</i> test results <sup>a</sup>
<i>L. monocytogenes</i>	1/2a	N14-2420	+
<i>L. monocytogenes</i>	1/2b	N14-2234	+
<i>L. monocytogenes</i>	1/2c	N14-2232	+
<i>L. monocytogenes</i>	3a	C23 FE8-1	+
<i>L. monocytogenes</i>	3c	N14-0326	+
<i>L. monocytogenes</i>	4b	N14-2079	+
<i>L. monocytogenes</i>	4-nm <sup>b</sup>	N14-0600	+
<i>L. aquatica</i>		DSM 26686	-
<i>L. cornellensis</i>		DSM 26689	-
<i>L. fleischmannii</i>		DSM 24998	-
<i>L. floridensis</i>		DSM 26687	-
<i>L. grandensis</i>		DSM 26688	-
<i>L. grayi</i>		DSM 20601	-
<i>L. innocua</i>		DSM 20649	-
<i>L. ivanovii</i>		DSM 20750	-
<i>L. riparia</i>		DSM 26685	-
<i>L. rocourtiae</i>		DSM 22097	-
<i>L. seeligeri</i>		DSM 20751	-
<i>L. weihenstephanensis</i>		DSM 24698	-
<i>L. welshimeri</i>		ATCC 35897	-

<sup>a</sup>Each isolate was incubated for 16 h in brain heart infusion (BHI) broth at 37°C; + or – indicates a positive or negative test result; <sup>b</sup>O-group 4 non-motile.

**Table 1:** Specificity of the Assurance GDS assays for *Listeria monocytogenes* and *Listeria* species.

### Detection limit of the Assurance GDS® *L. monocytogenes* assay

For these experiments, stationary phase BHI broth cultures of the seven *L. monocytogenes* strains (prepared as outlined above; about 10<sup>9</sup> cfu/ml) were 10-fold serially diluted in saline solution (0.85%) to obtain 10 ml broth cultures containing approximately 10<sup>3</sup>, 10<sup>4</sup>, 10<sup>5</sup>, 10<sup>6</sup>, 10<sup>7</sup>, and 10<sup>8</sup> cfu/ml. Subsets of each concentration and strain were then tested in triplicate using the Assurance GDS® *L. monocytogenes*

assay (Bio Control Systems) according to the manufacturers' instructions.

### Performance of the assurance GDS® *L. monocytogenes* assay in spiked food samples using a one-broth enrichment strategy

In the challenge test experiments, three ready-to-eat (RTE) food products were tested: steak tartare (dish from finely chopped or minced raw beef), bologna type sausage (traditional Swiss cooked sausage), and Gorgonzola cheese (veined Italian blue cheese). Tested food products were obtained from commercial retailers in Switzerland. In the laboratory, food samples were spiked with strain N14-2420 (*L. monocytogenes* serotype 1/2a). This strain was selected due to the frequent occurrence of serotype 1/2a strains in foods and their increasing proportion among human infections [14-18]. For each product type, three different spiking levels were used (10<sup>4</sup>, 10<sup>3</sup>, and <10<sup>2</sup> cfu/g). Spiked food samples (two for each spiking level and product type) were then enriched (25 g in 225 ml of Half-Fraser broth; Oxoid) for up to 48 h at 30°C. After 16 h of incubation in Half-Fraser broth (Oxoid), the Assurance GDS® *L. monocytogenes* assay (Bio Control Systems) was performed according to the manufacturers' instructions. In the case of a negative result, the Assurance GDS® test was repeated after 24 h and 48 h of incubation. *L. monocytogenes* colony counts of the spiked food samples and the enrichment broths (after 16 h, after 24 h, and if necessary after 48 h) were determined using Chromogenic *Listeria* agar (48 h at 37°C) supplemented with *Listeria* selective and differential supplement (Oxoid).

### Results and Discussion

The Assurance GDS® *L. monocytogenes* assay reliably detected all target strains in pure culture experiments, whereas isolates of other *Listeria* species yielded negative results (Table 1). Bosilevac et al. [12] reported for the first time specificity results for the Assurance GDS® test kit. However, in the present study, additional *L. monocytogenes* serotypes (3a and 3c) and a variety of recently described *Listeria* species (e.g. *L. aquatica*, *L. cornellensis*, *L. fleischmannii*, *L. floridensis*, *L. grandensis*, *L. riparia*, *L. rocourtiae*, *L. weihenstephanensis*) were also included.

To determine the detection limit of the Assurance GDS® *L. monocytogenes* assay, 10-fold serial dilutions of the seven *L. monocytogenes* strains (pure cultures; concentrations: 10<sup>3</sup>-10<sup>9</sup> cfu/ml) were tested in triplicate. Detailed evaluation results are shown in Table 2. With the exception of one run of the serotype 3c strain, the Assurance GDS® assay constantly detected the *L. monocytogenes* strains at concentrations ≥ 10<sup>6</sup> cfu/ml. A more heterogeneous picture was evident at 10<sup>5</sup> cfu/ml and 10<sup>4</sup> cfu/ml. Negative Assurance GDS® results in at least one run were observed at 10<sup>5</sup> cfu/ml for four strains and at 10<sup>4</sup> cfu/ml for six strains (four of them negative in two or three runs). At 10<sup>3</sup> cfu/ml, the Assurance GDS® assay did not detect any of the *L. monocytogenes* strains. Hence, concentrations >10<sup>5</sup> cfu/ml were required for reliable detection of *L. monocytogenes* using the Assurance GDS® system. This emphasizes the need for an adequate enrichment step (as specified by the manufacturer) when examining food products using this system.

For the challenge test experiments (steak tartare, bologna type sausage, Gorgonzola cheese; two different products of each type), RTE food samples were first spiked with the serotype 1/2a *L.*

*monocytogenes* strain (three different initial spiking levels) and then enriched in Half-Fraser broth (for up to 48 h).

	Serotype	Assurance GDS® <i>L. monocytogenes</i> test results at different concentrations (cfu/ml)						
		10 <sup>9</sup>	10 <sup>8</sup>	10 <sup>7</sup>	10 <sup>6</sup>	10 <sup>5</sup>	10 <sup>4</sup>	10 <sup>3</sup>
N14-2420	1/2a	+/+/+	+/+/+	+/+/+	+/+/+	+/+/-	+/-/-	-/-/-
N14-2234	1/2b	+/+/+	+/+/+	+/+/+	+/+/+	+/+/-	+/+/-	-/-/-
N14-2232	1/2c	+/+/+	+/+/+	+/+/+	+/+/+	+/-/-	-/-/-	-/-/-
C23 FE8-1	3a	+/+/+	+/+/+	+/+/+	+/+/+	+/+/+	+/+/+	-/-/-
N14-0326	3c	+/+/+	+/+/+	+/+/+	+/+/-	+/-/-	-/-/-	-/-/-
N14-2079	4b	+/+/+	+/+/+	+/+/+	+/+/+	+/+/+	+/+/-	-/-/-
N14-0600	4-nmb	+/+/+	+/+/+	+/+/+	+/+/+	+/+/+	+/+/-	-/-/-

<sup>a</sup>Each *L. monocytogenes* isolate was incubated for 16 h in brain heart infusion (BHI) broth at 37°C. After decimal serial dilution (10<sup>9</sup> to 10<sup>3</sup> cfu/ml), subsets of each concentration were tested in triplicate in the Assurance GDS® *L. monocytogenes* assay; each + or – indicates a positive or negative test result; <sup>b</sup>O-group 4 non-motile.

**Table 2:** Detection of *L. monocytogenes* by the Assurance GDS® *L. monocytogenes* assay at different concentrations (cfu/ml).

Product/spiking level	Assurance GDS® results and colony counts (cfu/ml) after enrichment <sup>a</sup>					
	After 16 h		After 24 h		After 48 h	
	GD S	cfu/ml	GD S	cfu/ml	GD S	cfu/ml
Tartare <sup>b</sup>						
<10 <sup>2</sup> cfu/g	–	1.9 x 10 <sup>2</sup>	+	2.7 x 10 <sup>4</sup>	nd	nd
	–	1.9 x 10 <sup>2</sup>	–	1.5 x 10 <sup>4</sup>	+	3.0 x 10 <sup>5</sup>
10 <sup>3</sup> cfu/g	+	2.2 x 10 <sup>7</sup>	nd	nd	nd	nd
	+	1.4 x 10 <sup>7</sup>	nd	nd	nd	nd
10 <sup>4</sup> cfu/g	+	8.8 x 10 <sup>7</sup>	nd	nd	nd	nd
	+	9.8 x 10 <sup>7</sup>	nd	nd	nd	nd
Bologna type sausage <sup>b</sup>						
<10 <sup>2</sup> cfu/g	–	2.0 x 10 <sup>2</sup>	+	3.0 x 10 <sup>3</sup>	nd	nd
	–	9.0 x 10 <sup>1</sup>	+	2.0 x 10 <sup>3</sup>	nd	nd
10 <sup>3</sup> cfu/g	+	3.0 x 10 <sup>6</sup>	nd	nd	nd	nd
	+	7.4 x 10 <sup>5</sup>	nd	nd	nd	nd
10 <sup>4</sup> cfu/g	+	4.2 x 10 <sup>7</sup>	nd	nd	nd	nd
	+	1.7 x 10 <sup>7</sup>	nd	nd	nd	nd
Gorgonzola <sup>b</sup>						
<10 <sup>2</sup> cfu/g	–	<10 <sup>2</sup>	na	2.0 x 10 <sup>2</sup>	+	4.2 x 10 <sup>4</sup>
	–	9.0 x 10 <sup>2</sup>	+	1.5 x 10 <sup>3</sup>	nd	nd
10 <sup>3</sup> cfu/g	+	1.0 x 10 <sup>6</sup>	nd	nd	nd	nd
	+	5.4 x 10 <sup>5</sup>	nd	nd	nd	nd

10 <sup>4</sup> cfu/g	+	4.0 x 10 <sup>6</sup>	nd	nd	nd	nd
	+	2.0 x 10 <sup>6</sup>	nd	nd	nd	nd

<sup>a</sup>Enrichment in Half-Fraser broth for up to 48 h at 30°C; + or – indicates a positive or negative test result; nd: not determined; na: no amplification; <sup>b</sup>Two different products were tested at each inoculation level.

**Table 3:** Results of the Assurance GDS® assay and colony counts after enrichment of food samples inoculated with *L. monocytogenes* (strain N14-2420, serotype 1/2a) in Half-Fraser broth.

Other studies addressing the performance of the Assurance GDS® test system for detection of *L. monocytogenes* in food products are so far widely lacking. Kerr and Bright [19] evaluated the Assurance GDS® assays for detection of *L. monocytogenes* and *Listeria* spp. in fish and seafood products by testing spiked and non-spiked samples after enrichment (for 18-22 h). These authors showed that the performance of the Assurance GDS® test systems were equivalent to the reference culture method, while being a much faster option.

In the present study, the Assurance GDS® assay reliably detected *L. monocytogenes* after 16 h of enrichment when the initial spiking level of the RTE food samples (steak tartare, bologna type sausage, Gorgonzola cheese) was between 10<sup>2</sup> and 10<sup>4</sup> cfu/g. The corresponding *L. monocytogenes* colony counts after 16 h of enrichment ranged from 5.4 x 10<sup>5</sup> to 9.8 x 10<sup>7</sup> cfu/ml in the enrichment broth (Table 3). On the other hand, the Assurance GDS® assay for *L. monocytogenes* yielded consistently negative results after 16 h of enrichment when the initial spiking levels were <10<sup>2</sup> cfu/g and the corresponding colony counts in the enrichment broth after 16 h <10<sup>3</sup> cfu/ml (Table 3). In these cases, incubation times of 24 h or even 48 h were required to obtain a positive Assurance GDS® test result. The corresponding *L. monocytogenes* colony counts (24 h or 48 h of enrichment and a positive Assurance GDS® test result) ranged from 1.5 x 10<sup>3</sup> to 3.0 x 10<sup>5</sup> cfu/ml. Thereby it must be considered that the *L. monocytogenes* counts on naturally contaminated (ready-to-eat) food products might vary widely, but frequently low contamination levels (<10<sup>2</sup> cfu/g) are expected [2,20-23].



Interestingly, one steak tartare sample showing  $1.5 \times 10^4$  cfu/ml after 24 h of enrichment was still negative in the GDS® Assurance *L. monocytogenes* assay, whereas enriched bologna type sausage and Gorgonzola cheese samples yielded positive test results when colony counts in the enrichment broth were in the range of  $1.5 \times 10^3$  to  $4.2 \times 10^4$  cfu/ml (Table 3). A special challenge was thereby the examination of Gorgonzola cheese samples with initial spiking levels  $<10^2$  cfu/g. The first two examined Gorgonzola samples yielded negative Assurance GDS® test results for *L. monocytogenes* even after 48 h of enrichment. Determination of colony counts showed the presence of a dominant *Listeria* spp. background flora, which hampered the growth of the spiked *L. monocytogenes*. Examinations were therefore repeated with two additional Gorgonzola samples that gave the results shown in Table 3.

## Conclusions

In summary, the Assurance GDS® *Listeria monocytogenes* assay has proven to be a reliable and easy to handle, rapid test system for the specific detection of *L. monocytogenes*. This system is suited as a tool for generating microbiological results, which can be used for a “positive batch release” especially also for RTE foods with short shelf lives. However, in a one-broth enrichment strategy, depending of the food matrix, enrichment times of 24 h or 48 h are required when the initial contamination level of the food matrix is  $<10^2$  cfu/g.

## References

- Allerberger F, Wagner M (2010) Listeriosis: a resurgent foodborne infection. Clin Microbiol Infect 16: 16-23.
- EFSA/ECDC (2015) The European Union summary report on trends and sources of zoonoses, zoonotic agents and food-borne outbreaks in 2014. EFSA Journal 13: 4329.
- Cartwright EJ, Jackson KA, Johnson SD, Graves LM, Silk BJ, et al. (2013) Listeriosis outbreaks and associated food vehicles, United States, 1998–2008. Emerg Infect Dis 19: 1-9.
- Ivanek R, Gröhn YT, Wiedmann M (2006) *Listeria monocytogenes* in multiple habitats and host populations: review of available data for mathematical modeling. Foodborne Pathog Dis 3: 319-336.
- Stephan R, Althaus D, Kiefer S, Lehner A, Hatz C, et al. (2015) Foodborne transmission of *Listeria monocytogenes* via ready-to-eat salad: a nationwide outbreak in Switzerland, 2013-2014. Food Control 57: 14-17.
- Swaminathan B, Gerner-Smidt P (2007) The epidemiology of human listeriosis. Microbes Infect 9: 1236-1243.
- Blatter S, Giezendanner N, Stephan R, Zweifel C (2010) Phenotypic and molecular typing of *Listeria monocytogenes* isolated from the processing environment and products of a sandwich-producing plant. Food Control 21: 1519-1523.
- Carpentier B, Cerf O (2011) Persistence of *Listeria monocytogenes* in food industry equipment and premises. Int J Food Microbiol 145: 1-8.
- Larivière-Gauthier G, Letellier A, Kérouanton A, Bekal S, Quessy S, et al. (2014) Analysis of *Listeria monocytogenes* strain distribution in a pork slaughter and cutting plant in the province of Quebec. J Food Prot 77: 2121-2128.
- ISO/EN 11290-1:2005 (2005) Microbiology of food and animal feeding stuffs - horizontal method for the detection and enumeration of *Listeria monocytogenes* - Part 1: Detection method (ISO 11290-1:1996 - AMD 1:2004).
- Jadhav S, Bhawe M, Palombo EA (2012) Methods used for the detection and subtyping of *Listeria monocytogenes*. J Microbiol Methods 88: 327-341.
- Bosilevac JM, Guerini MN, Koohmaraie M (2009) Increased detection of *Listeria* species and *Listeria monocytogenes* in raw beef, using the Assurance GDS molecular detection system with culture isolation. J Food Prot 72: 674-679.
- Brandão D, Liébana S, Pividori MI (2015) Multiplexed detection of foodborne pathogens based on magnetic particles. N Biotechnol 32: 511-520.
- Althaus D, Lehner A, Brisse S, Maury M, Tasara T, et al. (2014) Characterization of *Listeria monocytogenes* strains isolated during 2011–2013 from human infections in Switzerland. Foodborne Pathog Dis 11: 753-758.
- Ebner R, Stephan R, Althaus D, Brisse S, Maury M, et al. (2015) Phenotypic and genotypic characteristics of *Listeria monocytogenes* strains isolated during 2011–2014 from different food matrices in Switzerland. Food Control 57: 321-326.
- Gianfranceschi MV, D'Ottavio MC, Gattuso A, Bella A, Aureli P (2009) Distribution of serotypes and pulsotypes of *Listeria monocytogenes* from human, food and environmental isolates (Italy 2002–2005). Food Microbiol 26: 520-526.
- Lopez-Valladares G, Tham W, Parihar VS, Helmersson S, Andersson B, et al. (2014) Human isolates of *Listeria monocytogenes* in Sweden during half a century (1958–2010). Epidemiol Infect 142: 2251-2260.
- Martín B, Perich A, Gómez D, Yangüela J, Rodríguez A, et al. (2014) Diversity and distribution of *Listeria monocytogenes* in meat processing plants. Food Microbiol 44: 119-127.
- Kerr D, Bright G (2014) Comparative validation study to demonstrate the detection of *Listeria* spp. and *Listeria monocytogenes* in fish and seafood products with Assurance GDS® for *Listeria* spp. and *Listeria monocytogenes*. International Association for Food Protection (IAFP) Annual Meeting 2014, Indianapolis.
- EFSA (2013) Analysis of the baseline survey on the prevalence of *Listeria monocytogenes* in certain ready-to-eat foods in the EU, 2010–2011. Part A: *Listeria monocytogenes* prevalence estimates. EFSA Journal 11 (6): 3241.
- Lianou A, Sofos JN (2007) A review of the incidence and transmission of *Listeria monocytogenes* in ready-to-eat products in retail and food service environments. J Food Prot 70: 2172-2198.
- Little CL, Sagoo SK, Gillespie IA, Grant K, McLauchlin J (2009) Prevalence and level of *Listeria monocytogenes* and other *Listeria* species in selected retail ready-to-eat foods in the United Kingdom. J Food Prot 72: 1869-1877.
- Thisted Lambertz ST, Nilsson C, Brådenmark A, Sylvén S, Johansson A, et al. (2012) Prevalence and level of *Listeria monocytogenes* in ready-to-eat foods in Sweden 2010. Int J Food Microbiol 160: 24-31.