

Departement für Kleintiere  
Klinik für Kleintiermedizin  
der Vetsuisse-Fakultät Universität Zürich

Direktorin: Prof. Dr. med. vet. Claudia Reusch, Dipl. ECVIM-CA

Arbeit unter Leitung von PD Dr. med. vet. Felicitas Boretti, Dipl. ACVIM & ECVIM-CA

**Characterization of recombinant human and bovine thyroid-stimulating hormone preparations by mass spectrometry and determination of their endotoxin content**

Inaugural-Dissertation

zur Erlangung der Doktorwürde der  
Vetsuisse-Fakultät Universität Zürich

vorgelegt von

**Sandra Schäfer**

Tierärztin  
von Lupsingen BL, Schweiz

genehmigt auf Antrag von

PD Dr. med. vet. Felicitas Boretti, Dipl. ACVIM & ECVIM-CA, Hauptreferentin

Zürich 2014

## **Inhaltsverzeichnis**

Abstract	3
Zusammenfassung	4
Manuskript	5 - 19
Additional 1-3	
Lebenslauf	

Characterization of recombinant human and bovine thyroid-stimulating hormone preparations  
by mass spectrometry and determination of their endotoxin content

The thyroid-stimulating hormone (TSH) stimulation test to confirm canine hypothyroidism is commonly performed using a recombinant human TSH (rhTSH), as up to date, canine TSH is not yet commercially available. Limiting factors for the use of rhTSH are its high costs and occasional difficulties in product availability. Less expensive bovine TSH preparations (bTSH) purified from bovine pituitary glands are readily commercially available. The aim of this study was to evaluate two different bTSH products as alternative to rhTSH using mass spectrometry.

More than 50 proteins, including other pituitary hormones, bovine albumin, hemoglobin, and tissue proteins were identified in the bTSH preparations. In contrast, rhTSH proved to be a highly pure product. Significantly higher endotoxin levels could be detected in all bTSH products compared to the rhTSH.

Both bTSH products are crude mixtures and therefore not an acceptable alternative to rhTSH. Their use should be discouraged to prevent unintended side effects.

Keywords: Bovine TSH, Recombinant human TSH, Mass spectrometry, Endotoxin

Characterization of recombinant human and bovine thyroid-stimulating hormone preparations  
by mass spectrometry and determination of their endotoxin content

Zur Diagnosesicherung einer Hypothyreose beim Hund wird üblicherweise ein TSH-Stimulationstest durchgeführt. Hierzu wird heute rekombinantes humanes TSH (rhTSH) verwendet, da bis dato kanines TSH kommerziell nicht erhältlich ist. Einschränkungen beim Gebrauch von rhTSH sind die hohen Kosten und gelegentliche Schwierigkeiten beim Erwerb für den veterinärmedizinischen Gebrauch. Bovine TSH (bTSH) Präparate, welche aus bovinen Hypophysen gewonnen und gereinigt werden, sind demgegenüber einfach und kostengünstig erhältlich. Das Ziel der vorliegenden Studie war die massenspektrometrische Untersuchung zweier verschiedener bTSH Produkte und der Vergleich mit rhTSH, um deren Einsatz als Alternative zum rhTSH zu prüfen. Ausserdem sollte bei allen drei Produkten mittels eines Limulus-Amöbozyten-Lysat Tests der Endotoxin-Gehalt untersucht werden.

In beiden bTSH Präparaten konnten mehr als 50 verschiedene Proteine identifiziert werden, darunter andere hypophysäre Hormone, bovines Albumin, Hämoglobin und einige Gewebeproteine. Demgegenüber war das rhTSH ein hoch reines Produkt. Im Weiteren wurde in den bTSH Präparaten ein signifikant höherer Endotoxingehalt gemessen als im rhTSH.

Unsere Untersuchungen belegen, dass beide bTSH Präparate unreine Stoffgemische und somit keine akzeptable Alternative zum rhTSH darstellen. Um ungewollte Nebenwirkungen zu vermeiden, muss von deren Verwendung abgeraten werden.

Stichwörter: Bovines TSH, Rekombinantes humanes TSH, Massenspektrometrie, Endotoxin

# **Characterization of recombinant human and bovine thyroid-stimulating hormone preparations by mass spectrometry and determination of their endotoxin content**

Sandra Schaefer<sup>1</sup>  
Email: [sschaefer@vetclinics.uzh.ch](mailto:sschaefer@vetclinics.uzh.ch)

Paul O Hassa<sup>2</sup>  
Email: [hassa@vetbio.uzh.ch](mailto:hassa@vetbio.uzh.ch)

Nadja S Sieber-Ruckstuhl<sup>1</sup>  
Email: [nsieber@vetclinics.uzh.ch](mailto:nsieber@vetclinics.uzh.ch)

Marion Piechotta<sup>3</sup>  
Email: [Marion.Piechotta@tiho-hannover.de](mailto:Marion.Piechotta@tiho-hannover.de)

Claudia E Reusch<sup>1</sup>  
Email: [creusch@vetclinics.uzh.ch](mailto:creusch@vetclinics.uzh.ch)

Bernd Roschitzki<sup>4</sup>  
Email: [bernd.roschitzki@fgcz.uzh.ch](mailto:bernd.roschitzki@fgcz.uzh.ch)

Felicitas S Boretti<sup>1\*</sup>  
\* Corresponding author  
Email: [fboretti@vetclinics.uzh.ch](mailto:fboretti@vetclinics.uzh.ch)

<sup>1</sup> Clinic for Small Animal Internal Vetsuisse Faculty, University of Zurich, Zurich, Switzerland

<sup>2</sup> Institute of Veterinary Biochemistry and Molecular Biology, Vetsuisse Faculty, University of Zurich, Zurich, Switzerland

<sup>3</sup> Clinic for Cattle, Endocrinology Laboratory, University of Veterinary Medicine Hannover, Hanover, Germany

<sup>4</sup> Functional Genomics Center Zurich (FGCZ), University of Zurich/ETH Zurich, Zurich, Switzerland

# **Abstract**

## **Background**

The TSH stimulation test to confirm canine hypothyroidism is commonly performed using a recombinant human TSH (rhTSH), as up to date, canine TSH is not yet commercially available. Limiting factors for the use of rhTSH are its high costs and occasional difficulties in product availability. Less expensive bovine TSH preparations (bTSH) purified from bovine pituitary glands are readily commercially available. The aim of this study was to evaluate two different bTSH products as alternative to rhTSH using mass spectrometry.

## **Results**

More than 50 proteins, including other pituitary hormones, bovine albumin, hemoglobin, and tissue proteins were identified in the bTSH preparations. In contrast, rhTSH proved to be a highly pure product. Significantly higher endotoxin levels could be detected in all bTSH products compared to the rhTSH.

## **Conclusions**

Both bTSH products are crude mixtures and therefore not an acceptable alternative to rhTSH. Their use should be discouraged to prevent unintended side effects.

## **Keywords**

Bovine TSH, Recombinant human TSH, Mass spectrometry, Endotoxin

## **Background**

The TSH stimulation test has long been recognized as an accurate examination of thyroid function for the diagnosis of hypothyroidism in dogs. Increases in T<sub>4</sub> (thyroxin) after the administration of exogenous TSH provide an assessment of the functional reserve capacity of the thyroid gland, and this helps discriminate true hypothyroidism from other conditions with low T<sub>4</sub> secretion [1,2]. Canine TSH would be the ideal substance to perform the test, and although canine recombinant TSH could successfully be synthesized, it is not yet commercially available [3,4]. Therefore, the test is commonly performed using a recombinant human TSH (rhTSH) preparation [5-8]. Limiting factors, however, are the high cost of rhTSH and the occasional difficulties in obtaining the product. TSH purified from bovine pituitary glands (bTSH), which is easily commercially available, was widely used before rhTSH had been introduced in veterinary medicine. Equivalent biological activity of a bTSH product and the rhTSH has been shown in healthy beagle dogs [9]. However, in veterinary but also human medicine side effects after the use of bTSH were occasionally observed including allergic and anaphylactoid reactions [10-14]. Although they were assumed to result from impurities contained in the bTSH, a detailed characterization of commercially available bTSH to confirm this assumption has never been performed.

Mass spectrometry (MS) has been used to monitor protein purifications and to identify components in protein or peptide mixtures [15-18]. MS has also been suggested as an excellent analytical tool to establish key differences between recombinant proteins and their

natural counterparts, especially when these products are intended for clinical use [17]. To our knowledge, no MS studies have been conducted to describe the composition of commercial bTSH products.

Important contaminants found in parenteral drugs are endotoxins, also known as lipopolysaccharides (LPSs), components of the outer membrane of gram-negative bacteria. LPS in the blood circulation can lead to systemic inflammation and endotoxin shock [19]. Possible sources of endotoxins in drugs are chemicals, raw materials, or equipment used in the preparation and purification of the products. The Limulus Amebocyte Lysate (LAL) test is an officially accepted and very sensitive method to test endotoxin levels. To date, there are no published studies that compare the endotoxin content of different TSH preparations.

Knowledge of the composition of bTSH and the identification of potentially dangerous contaminants, e.g., endotoxins in this product are both prerequisites for veterinary clinical use of bTSH.

Therefore, aim of this study was to characterize commercially available purified bTSH (two different products) and rhTSH (one product) by MS. In addition, the levels of endotoxin were measured in both preparations by the Limulus endotoxin assay.

## Results

### Peptide and protein identification in bTSH

bTSH and rhTSH samples were separated by SDS-PAGE before analyses by MS. Selected examples of Bis-Tris gels after electrophoresis and Coomassie Blue staining of bTSH and rhTSH-loaded gels are shown in Figure 1A (bTSH product 1, Sigma Aldrich) and 1B (rhTSH), respectively.

---

**Figure 1 Bis-Tris SDS PAGE gel of bTSH (product 1) and rhTSH (each 70 µg/lane); after electrophoretic separation, protein bands were detected by Coomassie Blue staining.** The right lane of each gel represents the molecular weight marker (Prestained Protein Ladder; Page Ruler, Fermentas; Thermo Scientific, Glen Burnie, MD, USA), the left lane of gel **A** and **B**, the bTSH and rhTSH, respectively.

---

In addition to bTSH, more than 50 other proteins were identified in all bTSH samples, including other pituitary hormones such as bovine prolactin, LH, arginine vasopressin (AVP), proopiomelanocortin (POMC), GH, oxytocin-neurophysin, and glycoprotein hormone alpha chain. Moreover, bovine albumin, hemoglobin, and cytoskeleton proteins were found in many of the tested samples. Interestingly, major differences in the composition of albumin were observed among the different lots of bTSH (bTSH product 1, Sigma Aldrich), while the pituitary hormones (prolactin, LH, AVP, glycoprotein hormone alpha chain, and TSH) were not significantly different between the lots. A list of selected proteins of both preparations is presented in Table 1 and 2, respectively. A complete summary of all identified proteins in bTSH product 1 (three different lots) and 2 (one lot) are presented as Additional file 1: Table S1 and Additional file 2: Table S2, respectively; accession number and number of assigned spectra are given.

**Table 1 Selected proteins identified in the bTSH product 1 (Thyrotropic hormone from bovine pituitary, Sigma Aldrich; 3 lot numbers) by database search following mass spectrometry**

Identified Proteins	LOT 069 K1588	LOT 119 K1583	LOT 040 M1246
Prolactin	143	105	60
Lutropin subunit beta	100	84	86
Vasopressin-neurophysin2-copeptin	85	87	61
Thyrotropin subunit beta	76	64	44
Glycoprotein hormones alpha chain	57	73	75
Pro-opiomelanocortin	36	35	42
Serum albumin	25	1	0
Somatotropin	12	11	8
Beta-2-glycoprotein 1	5	2	0
Hemoglobin subunit beta	4	2	3
Oxytocin-neurophysin 1	4	2	2

Numbers of assigned spectra are given (corresponding to the numbers of the acquired spectra per protein).

**Table 2 Selected proteins identified in the bTSH product 2 (TSH, bovine pituitary, Calbiochem Merck; 1 lot number) by database search following mass spectrometry**

Identified Proteins	LOT
	D00106386
Beta-2-glycoprotein 1	52
Serum albumin	51
Lactotransferrin	45
Hemoglobin subunit beta	32
Lutropin subunit beta	16
Glycoprotein hormones alpha chain	13
Thyrotropin subunit beta	13



Vasopressin-neurophysin 2-copeptin	7
Somatotropin	6
Pro-opiomelanocortin	4
Hemoglobin subunit alpha	3
Follitropin subunit beta	3

Numbers of assigned spectra are given (corresponding to the numbers of the acquired spectra per protein).

## Peptide and protein identification in rhTSH

In the rhTSH only TSH and several enzymes (e.g., caspase, galectin, serpin, calmodulin-like protein, and glyceraldehyde-3-phosphate dehydrogenase species) could be identified. In contrast to bTSH, the composition of rhTSH was extremely consistent between the different lot numbers.

A list of selected proteins from the two lot numbers of rhTSH is shown in Table 3 and is given as the number of assigned spectra. A complete summary of all identified proteins in the rhTSH products is presented as Additional file 3: Table S3; accession number and number of assigned spectra are given.

**Table 3 Selected proteins identified in the rhTSH preparation (Thyrogen, Genzyme GmbH; 2 lot numbers) by database search following mass spectrometry**

Identified Proteins	LOT	LOT
	A8035H40	A8063H19
Thyrotropin subunit beta	116	126
Glycoprotein hormones alpha chain	25	26

Numbers of assigned spectra are given (corresponding to the numbers of the acquired spectra per protein).

## Estimation of pituitary protein level

The number of detectable peptides of the identified pituitary proteins was comparable. Therefore it was assumed that the number of the detected spectra reflected the individual abundances of the proteins, which allowed an estimation of the hormone quantities. Based on this assumption, the most commonly occurring proteins in the bTSH product 1 (Sigma Aldrich) were prolactin and LH and although not statistically significant there was a trend towards higher levels of both hormones compared to the TSH ( $p = 0.1$ ) and a trend towards lower GH levels compared to those of TSH, LH and prolactin (Figure 2). The rhTSH preparation had only minor amounts of proteins other than TSH.

**Figure 2 Numbers of assigned spectra of the hormones detected in the bTSH (product 1, Sigma Aldrich).** Each symbol represents a lot number. The line represents the median of the 3 different lot numbers of each hormone. Thyroid stimulating hormone (TSH), luteinizing hormone (LH), arginine vasopressin (AVP), proopiomelanocortin (POMC), growth hormone (GH).

Contaminants such as keratin or trypsin, and proteins that were detected in bTSH and in rhTSH were eliminated from the analyses.

## Measurement of LH and GH by ELISA

LH and GH were additionally measured by ELISA to confirm the findings of the MS with respect to the quantity of the hormone levels. The LH levels in the bTSH (product 1, Sigma Aldrich) samples were significantly higher compared to GH levels. Variations between the three different bTSH lots were only minimal and not statistically significant. Both, GH and LH levels in the rhTSH preparations were undetectable (Table 4).

**Table 4 Bovine luteinizing hormone (LH) and bovine growth hormone (GH) concentrations in bTSH and rhTSH preparations as determined by ELISA**

TSH product	LOT	LH (µg/ml)	GH (µg/ml)
bTSH	040 M1246	34.9	2.59
		37.4	2.29
	119 K1583	35.4	2.97
		34.5	2.75
rhTSH	A8063H19	< 1.2 ng/ml	< 0.8 ng/ml

## Endotoxin levels

Results of the LAL kinetic chromogenic assays of bTSH and rhTSH preparations at different dilutions are presented in Table 5. Endotoxin levels of all bTSH preparations were significantly higher compared to rhTSH. Further, large variations in endotoxin levels were observed among the different lots of the bTSH product 1 (Sigma Aldrich).

**Table 5 Endotoxin concentrations in two bTSH products (Sigma Aldrich, 3 lot numbers; Calbiochem Merck 1 lot number) and rhTSH (2 lot numbers) as determined by Limulus amoebocyte lysate (LAL) kinetic chromogenic assay**

Sample	LOT	Endotoxin/ml	Endotoxin/mg
bTSH (Sigma Aldrich)	069 K1588	434 EU/ml	80.4 EU/mg
	119 K1583	158.4 EU/ml	29.3 EU/mg
	040 M1246	133.6 EU/ml	24.7 EU/mg

<b>bTSH</b>			
	D00106386	622 EU/ml	12.4 EU/mg
(Calbiochem Merck)			
<b>rhTSH</b>	A8035H40	0.15 EU/ml	0.14 EU/mg
(Genzyme GmbH)	A8063H19	<0.05 EU/ml	<0.04 EU/mg

## Discussion

Mass spectrometry analyses demonstrated that all bTSH samples of two different products contained in addition to TSH, more than 50 other proteins, including not only other pituitary hormones (e.g. LH, prolactin, AVP, POMC, and GH) but also bovine albumin, hemoglobin, and several tissue proteins (e.g., collagen, superoxide dismutase, and cathepsin). In contrast, in the rhTSH samples only TSH and a low number of enzymes used for its production were detected.

Recombinantly produced proteins usually contain only the protein or peptide coded by the target gene, and some minor enzymes used in the commercial production process. rhTSH (Thyrogen) is produced in genetically engineered cell lines (Chinese hamster ovary cells-CHO) and is highly purified by a combination of ion exchange and dye affinity chromatography [20] (Genzyme Corporation, [www.genzyme.com](http://www.genzyme.com)).

In contrast to rhTSH, the bTSH products are purified proteins, prepared from the entire bovine pituitary gland. Depending on the purification procedure, it may be difficult to obtain a real pure product, particularly if proteins with high similarities to the target protein (in our case TSH) are present. It was therefore not surprising that other pituitary hormones were found in the various bTSH lots. Major cell types of the anterior pituitary are somatotrophic cells, constituting approximately 50% of the cell population, followed by lactotropic (10-25%), corticotrophic (10-20%), gonadotropic (10%), and thyrotrophic cells (10%) [21].

Based on the high proportion of somatotrophic cells in a pituitary mixture, a higher quantity of growth hormone would be expected if the protein purification method was non-specific. However, the proportion of GH in bTSH as determined by MS analyses and ELISA was low, confirming that the TSH purification was to some extent specific.

TSH belongs to the glycoprotein hormone family that includes FSH, hCG and LH; they are heterodimeric proteins and consist of a common  $\alpha$ -subunit (or  $\alpha$ -chain) and a unique  $\beta$ -subunit, which confers biological specificity to each hormone. Due to the very high homology between TSH and LH, separation of the two hormones is difficult and necessitates extensive purification procedures [22-24], which in turn is laborious and expensive. Based on our results showing a rather high proportion of LH (demonstrated by MS and ELISA), it can be assumed that such a special purification procedure has not been applied in the case of our bTSH.

In a previous study, we found that the administration of 75  $\mu$ g rhTSH to healthy beagle dogs resulted in stimulation of T4 similar to that caused by administration of 500  $\mu$ g bTSH (product1) [9]. Based on the current results, the more than six times higher dose of bTSH required to elicit the same effect as rhTSH can be explained by the high content of

contaminating proteins in the bTSH lots. As mentioned above, the amounts of GH in the bTSH products were low, which led us to assume that at least some specific purification procedure was used; the more surprising was the finding of the highly variable amounts of bovine albumin that we could detect. Before rhTSH was available, bTSH given to human patients caused allergic reactions in up to 43% of the cases; some patients required emergency treatment [12-14]. Contaminating proteins such as albumin or globulins and antibodies against these proteins were suspected to be a likely cause and most side effects were observed after repeated injections of bTSH [12,13]. In veterinary medicine, side effects after the use of bTSH include allergic and anaphylactoid reactions [10,11]. Although these reactions were assumed to result from impurities contained in the bTSH, a detailed characterization to confirm this assumption has not been performed. If foreign albumin is administered, species differences in albumin can lead to anaphylactic and immunologic reactions. This has been shown both in human patients and in dogs [25-30]. Clinical symptoms in dogs receiving non-canine albumin ranged from edema (facial, distal limbs), urticaria, vomiting, and diarrhea to severe shock-like reactions (hypotension, collapse). Some of these symptoms have also been observed in dogs receiving bTSH, making adverse reactions to albumin contained in bTSH a likely factor. In one veterinary study, reactions had only been observed after repeated administration of bTSH. The authors hypothesized that these reactions were the result of hypersensitization [11]. Clearly, to avoid hypersensitization, a canine TSH would be the most appropriate substance for performing the TSH stimulating test in dogs. However, to our knowledge, purified canine TSH is not commercially available in sufficient amounts and although recombinant canine TSH could successfully be synthesized, it is not yet on the market [3,4].

However, in one report, two fatal outcomes were described after a single administration of bTSH to dogs that had no prior exposure to bTSH [10]. A possible explanation in these cases, in which a hypersensitization can be excluded, would be the high endotoxin concentration in some bTSH lots. Based on the LAL assay, the endotoxin concentration in bTSH was up to 500 times higher than that measured in rhTSH (434 and 0.27 EU/mL, respectively). Endotoxins in the blood circulation can lead to systemic inflammation and endotoxin shock with severe cardiovascular disturbances, multi-organ failure, and possible fatal outcome [19]. The highly variable endotoxin content in different lots makes adverse reactions unpredictable, and hence, the administration of bTSH to dogs is risky. Based on the high level of purity and the lack or minimal amount of endotoxin content, adverse reactions to rhTSH are unlikely. Therefore, compared to bTSH, rhTSH can be regarded as a much safer product. This assumption is supported by the finding that thus far, no adverse reactions in dogs have been observed even after the repeated administration of rhTSH [5,6,9]. In selected dogs in which hypothyroidism had been suspected, we repeated the TSH stimulation test three times using rhTSH. Even though the third test was performed within 12 to 24 months after the first stimulation, no dogs had adverse reactions after the rhTSH administration [6].

The bTSH is commonly used also for *in vitro* experiments to test the influence of TSH [31-34]. With our data, we could show that other pituitary hormones, albumin, endotoxins and an extensive list of contaminating proteins (e.g., collagen, hemoglobin, superoxide dismutase, and cathepsin) are constituents in two different bTSH products. The effects of these impurities are rarely if ever considered in interpretations of the results from these cell culture experiments. However, in the light of our findings, the rhTSH should be preferred for *in vitro* usage if the sole effect of TSH has to be tested.

## Conclusion

The results of the present study show that rhTSH is a very pure product with minimal to no amounts of endotoxin, whereas bTSH is a crude mixture of proteins with a high level of endotoxin. Both, contaminating proteins like bovine albumin and the high endotoxin content in the bTSH can potentially lead to adverse reactions, which have been observed in the past, not only in humans, but also in veterinary medicine [10-14]. The highly variable and lot-dependent composition of the bTSH makes adverse reactions unpredictable. Consequently, although less expensive and easily available, bTSH cannot be regarded as an alternative to the more expensive rhTSH, and its use should be discouraged to prevent unintended side effects *in vitro* and *in vivo*.

## Methods

### Sample preparation

Two different bovine TSH products (Thyrotropic hormone from bovine pituitary; product 1 Sigma Aldrich, Buchs, Switzerland and product 2 Calbiochem, Merck Chemicals, Lucerne, Switzerland) and rhTSH (Thyrogen, Genzyme GmbH, Baar, Switzerland) were analyzed. All TSH products were dissolved in aqua ad iniectabilia (Fresenius Kabi AG, Stans, Switzerland), and 70 µg total protein of each TSH sample were loaded separately onto a 12%- Bis-Tris gel (NuPage Novex Bis-Tris Mini Gel, Invitrogen AG, Basel, Switzerland). After separation, the gels were stained with colloidal Coomassie Blue (RotiBlue, Roth, Karlsruhe, Germany) and the stained gel lanes were excised, digested in-gel, and finally dissolved in 3% acetonitrile/ 0.1% formic acid for MS analyses as previously described [35].

### Mass spectrometry analyses

Analyses were performed with an LTQ-FT-ICR Ultra mass spectrometer (Thermo Fisher Scientific, Pittsburgh, PA, USA) equipped with an Eksigent-Nano-high pressure liquid chromatography system (Eksigent Technologies, Dublin, CA, USA) as described by Moretti et al. [35]. An inclusion list was used containing the m/z values of all doubly and triply charged tryptic peptides from human (P01222) and bovine (P01223) TSH.

### Peptide and protein identification (data analyses)

The MS data were analyzed and peptides identified with Mascot (Version 2.3.0, Matrix Science) using the Swiss-Prot database (release date, 05.Oct10; 521016 sequence entries). The peptide tolerance was set to  $\pm 5$  ppm, the MS/MS tolerance to  $\pm 0.6$  Da, and carbamidomethylation of cysteine was specified as a fixed modification.

### Determination of bovine growth hormone concentration

The concentration of bovine growth hormone (GH) was measured in duplicate by an Enzyme-linked Immunosorbant Assay (ELISA) using a GH antibody (10 mg/mL, 100 µL total volume, 1: 200 000 dilution, anti-ovine GH-3, Lot# AFP0802210Rb) obtained from Dr. Parlow [National Hormone & Peptide Program (NHPP) of the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK), National Institutes of Health, Bethesda, MD, USA]. The assay was performed as previously described, but with the following

modifications [36,37]. The anti-ovine GH antibody was distributed into all wells of microtiter plates, which had been coated with anti-rabbit-globulin antiserum, and incubated for 24 h at room temperature. After decanting the supernatants, 15  $\mu$ L of different concentrations of GH standard (0.8–100 ng/mL, NIDDK-bGH, AFP-10325C) diluted in assay buffer (7.1 g/L  $\text{Na}_2\text{HPO}_4$ , 1.1 g/L  $\text{KH}_2\text{PO}_4$ , 1.2 g/L NaCl, and 1.8 g/L EDTA, pH 7.5; containing 1% chicken serum; 100  $\mu$ g/well), and samples were added and incubated for 24 h. After incubation, the supernatants were discarded and biotin-labeled GH was added and incubated for 3 h at room temperature. Finally, the substrate was added (40 min at 37°C), and after stopping, the reaction, the optical density was measured at 450 nm. Only concentrations in the linear range of the standard curve were evaluated. Therefore, lower and upper detection limit for undiluted samples were 2.0 and 50.0 ng/mL, respectively (working range). The EIA method was validated for protein extraction samples by analyzing linearity of dilution and intra-assay coefficient of variation (CV) using the respective TSH diluted in GH assay buffer. The signal linearity of dilution was 81–96%, the intra-assay CV was 15.1%, and the 50% binding of this assay system was 6.2 ng/mL. Protein samples were diluted 1:10 in water and additionally 1:20 in peptide buffer.

### **Determination of bovine luteinizing hormone concentration**

Bovine luteinizing hormone (LH) concentration was measured in duplicate using a solid phase two-site enzyme immunoassay kit (LH Detect, INRA, Tours, France). The immunoassay was validated for purified protein samples by analyzing linearity of dilution, and the intra-assay CV. The assay sensitivity was 1 ng/mL, and the intra-assay CV was 9.6%. Linearity of dilution was 82–102%. Protein for LH detection was diluted in water at 1:100 and additionally 1:100 in assay buffer. The standard curve ranged from 1.2 - 40 ng/mL.

### **Endotoxin assay**

Endotoxin levels were measured using a Limulus amoebocyte lysate (LAL) kinetic chromogenic assay (Pyrogen LAL Chromogenic Assay, Lonza, Verviers, Belgium). The test is based on the measurement of the chromophore released from a suitable chromogenic peptide by the reaction of endotoxins with the lysate. Bacterial endotoxin was measured with upper and lower detection limits of 0.05 and 50 Endotoxin units (EU)/mL, respectively. The sample was mixed with the LAL/substrate reagent, placed in an incubating plate reader, and monitored over time for the appearance of a yellow color, according to the instructions of the manufacturer. The same TSH lots were used for endotoxin assays and mass spectrometry. Results are given in both EU/mL and EU/mg.

### **Statistical analyses**

Statistical support was provided by the biostatistical staff of the FGCZ. MS data were analyzed using Scaffold 3 (Proteome Software, Portland, OR, USA). The threshold for positive identification of a given protein was set to at least two peptides with a peptide probability of 95% or better. The minimal sequence coverage was set to a minimum of 10%. Data were further analyzed using non-parametric statistical methods (SPSS, Statistical Package for the Social Science, Software Packets for Windows, Version 18 and GraphPad PRISM<sup>®</sup> for Windows, Version 5.0.) Kruskal-Wallis, Wilcoxon and Dunn's multiple comparisons test were used. The Mann-Whitney-U-Test was used to determine differences between two groups (hormones, endotoxins). Values of  $p < 0.05$  were considered statistically significant.

## Competing interests

None of the authors of this paper has a financial or personal relationship with other people or organizations that could inappropriately influence or bias the content of the paper. The authors declare that they have no competing interests.

## Authors' contributions

SS carried out all the experimental work of the MS and drafted the manuscript. MP carried out the immunoassays and helped with revision of the manuscript. POH and BR coordinated and supervised the experimental work and helped with revision of the manuscript. BR did the biostatistical analyses. NSSR and CER supervised and advised on scientific content of the manuscript and critical revision of the text. FB designed, coordinated and supervised the study and manuscript preparation. All authors read and approved the final manuscript.

## Acknowledgements

This work was supported by the Vetsuisse Faculty Research Program 'CLINOMICS', University of Zurich Switzerland and the Albert Heim Foundation for Canine Research, Bern, Switzerland. We thank Theres Meili and Anamarija Ivos for excellent technical assistance.

## References

1. Feldman EC, Nelson RW: **Hypothyroidism**. In *Canine and feline endocrinology and reproduction*. Edited by Feldman EC, Nelson RW. W.B.Saunders Company; 1996:68–117.
2. Ferguson DC: **Update on diagnosis of canine hypothyroidism**. *Vet Clin North Am Small Anim Pract* 1994, **24**(3):515–539.
3. Yang X, McGraw RA, Ferguson DC: **cDNA cloning of canine common alpha gene and its co-expression with canine thyrotropin beta gene in baculovirus expression system**. *Domest Anim Endocrinol* 2000, **18**(4):379–393.
4. Yanga X, McGraw RA, Su X, Katakam P, Grosse WM, Li OW, Ferguson DC: **Canine thyrotropin beta-subunit gene: cloning and expression in Escherichia coli, generation of monoclonal antibodies, and transient expression in the Chinese hamster ovary cells**. *Domest Anim Endocrinol* 2000, **18**(4):363–378.
5. Boretti FS, Sieber-Ruckstuhl NS, Favrot C, Lutz H, Hofmann-Lehmann R, Reusch CE: **Evaluation of recombinant human thyroid-stimulating hormone to test thyroid function in dogs suspected of having hypothyroidism**. *Am J Vet Res* 2006, **67**(12):2012–2016.
6. Boretti FS, Sieber-Ruckstuhl NS, Wenger-Riggenbach B, Gerber B, Lutz H, Hofmann-Lehmann R, Reusch CE: **Comparison of 2 doses of recombinant human thyrotropin for thyroid function testing in healthy and suspected hypothyroid dogs**. *J Vet Intern Med* 2009, **23**(4):856–861.


7. Daminet S, Fifle L, Paradis M, Duchateau L, Moreau M: **Use of recombinant human thyroid-stimulating hormone for thyrotropin stimulation test in healthy, hypothyroid and euthyroid sick dogs.** *Can Vet J* 2007, **48**(12):1273–1279.
8. Sauve F, Paradis M: **Use of recombinant human thyroid-stimulating hormone for thyrotropin stimulation test in euthyroid dogs.** *Can Vet J* 2000, **41**(3):215–219.
9. Boretti FS, Sieber-Ruckstuhl NS, Willi B, Lutz H, Hofmann-Lehmann R, Reusch CE: **Comparison of the biological activity of recombinant human thyroid-stimulating hormone with bovine thyroid-stimulating hormone and evaluation of recombinant human thyroid-stimulating hormone in healthy dogs of different breeds.** *Am J Vet Res* 2006, **67**(7):1169–1172.
10. Hasler A, Rohner K: **Serious reactions after the TSH-stimulation test in the dog.** *Schweiz Arch Tierheilkd* 1992, **134**(9):423–427.
11. Kobayashi DL, Nichols R, Peterson ME: **Serum thyroid hormone concentrations in clinically normal dogs after administration of freshly reconstituted vs previously frozen and stored thyrotropin.** *J Am Vet Med Assoc* 1990, **197**(5):597–599.
12. Kirkpatrick CH, Meek JC, Rich RR: **Mechanism of allergy to components of commercial bovine thyrotropin.** *J Allergy Clin Immunol* 1973, **51**(5):296–302.
13. Krishnamurthy GT: **Human reaction to bovine TSH: concise communication.** *J Nucl Med* 1978, **19**(3):284–286.
14. Sherman WB, Werner SC: **Generalized allergic reaction to bovine thyrotropin.** *JAMA* 1964, **190**:244–245.
15. Domon B, Aebersold R: **Mass spectrometry and protein analysis.** *Science* 2006, **312**(5771):212–217.
16. Domon B, Alving K, He T, Ryan TE, Patterson SD: **Enabling parallel protein analysis through mass spectrometry.** *Curr Opin Mol Ther* 2002, **4**(6):577–586.
17. Flensburg J, Belew M: **Characterization of recombinant human serum albumin using matrix-assisted laser desorption ionization time-of-flight mass spectrometry.** *J Chromatogr A* 2003, **1009**(1–2):111–117.
18. Jensen C, Haebel S, Andersen SO, Roepstorff P: **Towards monitoring of protein purification by matrix-assisted laser desorption ionization mass spectrometry.** *Int J Biochem Mass Spectrometry Ion Processes* 1997, **160**:339–356.
19. Galanos C, Freudenberg MA: **Mechanisms of endotoxin shock and endotoxin hypersensitivity.** *Immunobiology* 1993, **187**(3–5):346–356.
20. Braverman LE, Pratt BM, Ebner S, Longcope C: **Recombinant human thyrotropin stimulates thyroid function and radioactive iodine uptake in the rhesus monkey.** *J Clin Endocrinol Metab* 1992, **74**(5):1135–1139.




21. Yeung CM, Chan CB, Leung PS, Cheng CH: **Cells of the anterior pituitary.** *Int J Biochem Cell Biol* 2006, **38**(9):1441–1449.
22. Arora T, Vashistha N, Muralidhar K: **Chromatographic separation of thyroid stimulating hormone from luteinizing hormone.** *J Liq Chromatogr R T* 2010, **33**(5):680–692.
23. Hartree AS, Butt WR, Kirkham KE: **The separation and purification of human luteinizing and thyrotrophic hormones.** *J Endocrinol* 1964, **29**:61–69.
24. Stanton PG, Hearn MT: **High-performance liquid chromatography of amino acids, peptides and proteins. LXXV. Isolation of isohormones of bovine thyrotropin and lutropin.** *J Chromatogr* 1987, **397**:379–387.
25. Cohn LA, Kerl ME, Lenox CE, Livingston RS, Dodam JR: **Response of healthy dogs to infusions of human serum albumin.** *Am J Vet Res* 2007, **68**(6):657–663.
26. Francis AH, Martin LG, Haldorson GJ, Lahmers KK, Luther TY, Alperin DC, Hines SA: **Adverse reactions suggestive of type III hypersensitivity in six healthy dogs given human albumin.** *J Am Vet Med Assoc* 2007, **230**(6):873–879.
27. Mackensen A, Drager R, Schlesier M, Mertelsmann R, Lindemann A: **Presence of IgE antibodies to bovine serum albumin in a patient developing anaphylaxis after vaccination with human peptide-pulsed dendritic cells.** *Cancer Immunol Immunother* 2000, **49**(3):152–156.
28. Pagan JA, Postigo I, Rodriguez-Pacheco JR, Pena M, Guisantes JA, Martinez J: **Bovine serum albumin contained in culture medium used in artificial insemination is an important anaphylaxis risk factor.** *Fertil Steril* 2008, **90**(5):2013. e2017-2019.
29. Wehner-Caroli J, Schreiner T, Schippert W, Lischka G, Fierlbeck G, Rassner G: **Anaphylactic reaction to bovine serum albumin after embryo transfer.** *Fertil Steril* 1998, **70**(4):771–773.
30. Wuthrich B, Stern A, Johansson SG: **Severe anaphylactic reaction to bovine serum albumin at the first attempt of artificial insemination.** *Allergy* 1995, **50**(2):179–183.
31. Bell A, Gagnon A, Dods P, Papineau D, Tiberi M, Sorisky A: **TSH signaling and cell survival in 3T3-L1 preadipocytes.** *Am J Physiol Cell Physiol* 2002, **283**(4):C1056–C1064.
32. Bell A, Gagnon A, Sorisky A: **TSH stimulates IL-6 secretion from adipocytes in culture.** *Arterioscler Thromb Vasc Biol* 2003, **23**(12):e65–e66.
33. Curcio F, Ambesi-Impiombato FS, Perrella G, Coon HG: **Long-term culture and functional characterization of follicular cells from adult normal human thyroids.** *Proc Natl Acad Sci U S A* 1994, **91**(19):9004–9008.
34. Westermarck B, Karlsson FA, Walinder O: **Thyrotropin is not a growth factor for human thyroid cells in culture.** *Proc Natl Acad Sci U S A* 1979, **76**(4):2022–2026.

35. Moretti M, Grunau A, Minerdi D, Gehrig P, Roschitzki B, Eberl L, Garibaldi A, Gullino ML, Riedel K: **A proteomics approach to study synergistic and antagonistic interactions of the fungal-bacterial consortium *Fusarium oxysporum* wild-type MSA 35.** *Proteomics* 2010, **10**(18):3292–3320.
36. Kawashima C, Kida K, Hayashi KG, Amaya Montoya C, Kaneko E, Matsunaga N, Shimizu T, Matsui M, Miyake Y, Schams D, *et al*: **Changes in plasma metabolic hormone concentrations during the ovarian cycles of Japanese black and holstein cattle.** *J Reprod Dev* 2007, **53**(2):247–254.
37. Roh SG, Matsunaga N, Miyamoto A, Hidaka S, Hidari H: **Competitive enzyme immunoassay for bovine growth hormone.** *Endocr J* 1997, **44**(1):195–198

**Figure 1**



**Figure 2**



## Additional file 1

Table 1 (complete, including accession number)

Complete list of the proteins identified in the bTSH product 1 (Thyrotropic hormone from bovine pituitary, Sigma Aldrich; 3 lot numbers) by database search following mass spectrometry. Numbers of assigned spectra are given and the minimal sequence coverage was set to a minimum of 10%. Contaminations like keratin or trypsin, which were registered in bTSH as well as in rhTSH were excluded from analyses.

Identified Proteins	Accession Number	LOT 069K1588	LOT 119K1583	LOT 040M1246
Prolactin	sp P01239	143	105	60
Lutropin subunit beta	sp P04651	100	84	86
Vasopressin-neurophysin 2-copeptin	sp P01180	70	80	53
Thyrotropin subunit beta	sp P01223	70	54	44
Glycoprotein hormones alpha chain	sp P01217	57	73	75
Cathepsin D	sp P80209	46	38	48
Pro-opiomelanocortin	sp P01190	36	35	42
Stathmin	sp Q3T0C7	22	11	15
Serum albumin	sp P49822	17	0	0
Macrophage migration inhibitory factor	sp P80177	16	19	6
Somatotropin	sp P01246	12	11	8
Proteasome activator complex subunit 1	sp Q4U5R3	10	5	3
Prefoldin subunit 6	sp Q17Q89	9	0	6
Glucosylceramidase	sp Q2KHZ8	9	3	0
Cystatin-C	sp P01035	9	1	0
Selenium-binding protein 1	sp Q2KJ32	8	18	1
Vasopressin-neurophysin 2-copeptin	sp P01185	8	6	5
Cytochrome b-c1 complex subunit 7	sp P00129	8	4	3
Cathepsin H	sp Q3T0I2	8	2	3
Serum albumin	sp P02769	8	1	0
Ubiquitin	sp P62990	7	18	10
Prefoldin subunit 1	sp Q3SZE2	7	4	6
Ribonuclease UK114	sp Q3T114	7	4	4
Vasopressin-neurophysin 2-copeptin (Fragment)	sp P01181	7	1	3
Protein S100-A11	sp P24480	7	3	2
Cathepsin S	sp P25326	7	2	3
Lysozyme C, non-stomach isozyme	sp P80189	7	0	1
Thyrotropin subunit beta	sp P01222	6	10	0
Epididymal secretory protein E1	sp P79345	6	3	3
Fumarylacetoacetate hydrolase domain-containing protein 2	sp Q2KIB0	6	2	0
Phosphatidylethanolamine-binding protein 1	sp P13696	5	35	9

Connective tissue growth factor	sp O18739	5	1	5
Thioredoxin	sp O97680	5	1	0
Beta-2-glycoprotein 1	sp P17690	5	2	0
Palmitoyl-protein thioesterase 1	sp P45478	5	0	0
Acidic mammalian chitinase	sp Q95M17	5	0	0
Superoxide dismutase [Cu-Zn]	sp P00442	4	29	8
ATP synthase-coupling factor 6, mitochondrial	sp P13618	4	5	7
Hemoglobin subunit beta	sp P02070	4	2	3
Oxytocin-neurophysin 1	sp P01175	4	2	2
Cytochrome c oxidase subunit 5B, mitochondrial	sp P00428	4	2	0
Cytochrome c oxidase subunit 6B1	sp P00429	4	0	4
Insulin-like growth factor-binding protein 6	sp Q05718	4	0	2
Apolipoprotein A-I	sp P02648	4	0	0
Chromogranin-A	sp P05059	3	6	4
Mitochondrial import inner membrane translocase subunit Tim9	sp Q2KIV2	3	6	5
Proenkephalin-A	sp P01211	3	7	3
Coiled-coil domain-containing protein 58	sp A4FUI1	3	2	6
Mitochondrial import inner membrane translocase subunit Tim10	sp Q2NKR1	3	3	2
Beta-2-microglobulin	sp P01888	3	0	3
Spleen trypsin inhibitor I	sp P04815	3	0	2
Protein FAM136A	sp Q2HJI3	3	0	1
Transthyretin	sp O46375	3	0	0
Insulin-like growth factor-binding protein 7	sp Q16270	2	2	5
Glutathione peroxidase 1	sp P00435	2	0	7
Cathelicidin-4	sp P33046	2	2	1
Prefoldin subunit 5	sp Q8HYI9	2	2	1
COX assembly mitochondrial protein homolog	sp Q3SZM6	2	0	0
TP53-regulated inhibitor of apoptosis 1	sp O43715	1	4	2
Mitochondrial import inner membrane translocase subunit Tim13	sp Q9Y5L4	1	2	2
Acylphosphatase-1	sp P41500	1	6	0
Cytochrome c oxidase assembly protein COX19	sp A8E4L1	1	1	2
Coiled-coil-helix domain-containing protein 3, mitochondrial	sp Q5E9D3	0	2	3
Acyl-CoA-binding protein	sp P07107	0	2	0

## Additional file 2

Table 2 (complete, including accession number)

Complete list of identified proteins in the bTSH product 2 (TSH, bovine pituitary, Calbiochem Merck; 1 lot number) by database search following mass spectrometry. Numbers of assigned spectra are given and the minimal sequence coverage was set to a minimum of 10%. Contaminations like keratin or trypsin, which were registered in bTSH as well as in rhTSH were excluded from analyses.

Identified Proteins	Accession Number	LOT D00106386
Beta-2-glycoprotein 1	sp P17690	52
Serum albumin	sp P02769	51
Lactotransferrin	sp P24627	45
Annexin A1	sp P46193	36
Annexin A6	sp P79134	35
Annexin A2	sp A2SW69	33
Hemoglobin subunit beta	sp P02070	32
Filamin-A	sp P21333	32
Phospholipase D3	sp Q2KJJ8	29
Stress-induced-phosphoprotein 1	sp Q3ZBZ8	25
Peptidyl-prolyl cis-trans isomerase B	sp P80311	25
Annexin A3	sp Q3SWX7	25
Eukaryotic translation initiation factor 5A-1	sp P63241	22
Gelsolin	sp Q3SX14	20
Spectrin alpha chain, brain	sp P07751	20
Metalloproteinase inhibitor 2	sp Q9TRZ7	20
Hepatoma-derived growth factor	sp Q9XSK7	20
Spectrin beta chain, brain	sp Q01082	19
Thrombospondin-1	sp Q28178	19
Hemopexin	sp Q3SZV7	17
Prolargin	sp Q9GKN8	17
Lutropin subunit beta	sp P04651	16
Adenylate kinase 2, mitochondrial	sp P08166	16
Glycylpeptide N-tetradecanoyltransferase 1	sp O70310	14
Interleukin enhancer-binding factor 3	sp Q12906	14
Rho GDP-dissociation inhibitor 1	sp P19803	13
Glycoprotein hormones alpha chain	sp P01217	13
Adenylyl cyclase-associated protein 1	sp Q3SYV4	13
Calcyclin-binding protein	sp Q3T168	13
Thyrotropin subunit beta	sp P01223	13
GTP:AMP phosphotransferase, mitochondrial	sp P08760	13
Hepatoma-derived growth factor-related protein 3	sp Q923W4	13

Vinculin	sp P18206	12
Transgelin	sp Q9TS87	12
Gamma-glutamyl hydrolase	sp A7YWG4	11
Neurofascin	sp Q94856	11
Elongation factor 1-gamma	sp Q3SZV3	11
LDLR chaperone MESD	sp Q3T0U1	11
Glucosylceramidase	sp Q2KHZ8	10
Malate dehydrogenase, mitochondrial	sp Q32LG3	10
Peroxiredoxin-1	sp Q5E947	10
Vesicle-associated membrane protein-associated protein A	sp Q0VCY1	10
Neuroblast differentiation-associated protein AHNAK	sp Q09666	9
SUMO-conjugating enzyme UBC9	sp P63279	9
High mobility group protein B1	sp A9RA84	9
Phosphatidylinositol-binding clathrin assembly protein	sp Q55012	9
Phosphatidylethanolamine-binding protein 1	sp P13696	9
Aminoacyl tRNA synthase complex-interacting multifunctional protein 1	sp P31230	9
GDNF family receptor alpha-1	sp P56159	8
Transcription elongation factor A protein 1	sp Q29RL9	8
Microtubule-associated protein RP/EB family member 1	sp Q3ZBD9	8
Peptidyl-glycine alpha-amidating monooxygenase	sp P10731	8
DNA-(apurinic or apyrimidinic site) lyase	sp P23196	8
Cathepsin L1	sp P25975	8
Protein canopy homolog 3	sp Q0P5N1	8
Serotransferrin	sp Q29443	8
Fermitin family homolog 2	sp Q96AC1	8
Gamma-aminobutyric acid receptor-associated protein-like 2	sp P60519	8
Ribosome-recycling factor, mitochondrial	sp Q0VCQ4	8
Cystatin-C	sp P01035	7
Methyl-CpG-binding protein 2	sp P51608	7
Myosin light chain kinase, smooth muscle	sp Q28824	7
Vasopressin-neurophysin 2-copeptin	sp P01180	7
Insulin-like growth factor-binding protein 7	sp Q16270	7
Tropomodulin-2	sp Q9NZR1	7
Tetranectin	sp Q2KIS7	7
D-tyrosyl-tRNA(Tyr) deacylase 1	sp Q2T9V8	7
BTB/POZ domain-containing protein KCTD12	sp Q6WVG3	7
Nascent polypeptide-associated complex subunit alpha, muscle-specific form	sp P70670	6
U2 small nuclear ribonucleoprotein A'	sp P09661	6
Complement factor B	sp P81187	6
Laminin subunit alpha-4	sp Q16363	6
Cochlin	sp O43405	6
Thioredoxin	sp O97680	6
78 kDa glucose-regulated protein	sp P06761	6

Plasminogen	sp P06868	6
RNA-binding protein FUS	sp P35637	6
Lumican	sp Q05443	6
Partner of Y14 and mago	sp A6QPH1	6
Somatotropin	sp P01246	6
Metalloproteinase inhibitor 1	sp P20414	6
Factor XIIa inhibitor	sp P50448	6
Ribosome maturation protein SBDS	sp Q3SWZ6	6
Glucosamine 6-phosphate N-acetyltransferase	sp Q5RAL9	6
Peroxiredoxin-2	sp Q9BGI3	6
Complement C3	sp Q2UVX4	5
Microtubule-associated proteins 1A/1B light chain 3 beta 2	sp A6NCE7	5
Sorbin and SH3 domain-containing protein 2	sp O94875	5
GlutaminyI-tRNA synthetase	sp Q3MHH4	5
Alpha-enolase	sp Q9XSJ4	5
Non-muscle caldesmon	sp Q62736	5
Elongation factor 1-alpha 1	sp A2Q0Z0	5
Microtubule-associated protein RP/EB family member 3	sp Q5XIT1	5
Clathrin coat assembly protein AP180	sp O60641	5
Heterogeneous nuclear ribonucleoprotein C	sp O77768	5
Heparin-binding growth factor 2	sp P03969	5
Haptoglobin	sp Q2TBU0	5
Calcium-dependent secretion activator 1	sp Q9ULU8	5
Eukaryotic translation initiation factor 3 subunit G	sp O75821	5
Nucleolin	sp P19338	5
Protein disulfide-isomerase A3	sp P38657	5
Ubiquitin carboxyl-terminal hydrolase 14	sp P54578	5
Coiled-coil-helix-coiled-coil-helix domain-containing protein 7	sp Q17Q91	5
Protein NipSnap homolog 3A	sp Q5RAA9	5
Protein dpy-30 homolog	sp Q2NKU6	5
Leiomodin-1	sp Q8BVA4	5
Spectrin beta chain, brain 1	sp Q62261	5
Heat shock cognate 71 kDa protein	sp A2Q0Z1	4
Protein AMBP	sp P00978	4
Lamina-associated polypeptide 2, isoform alpha	sp P42166	4
Transcription factor BTF3	sp P20290	4
Collagen alpha-1(XIV) chain	sp Q80X19	4
Pro-opiomelanocortin	sp P01190	4
Proteasomal ubiquitin receptor ADRM1	sp A1L5A6	4
Transcription activator BRG1	sp A7Z019	4
Collagen alpha-1(II) chain	sp P02458	4
Glia-derived nexin	sp P07093	4
40S ribosomal protein S19	sp P17074	4
60S ribosomal protein L12	sp P23358	4



40S ribosomal protein S3	sp P23396	4
Transgelin-2	sp P37802	4
Mesencephalic astrocyte-derived neurotrophic factor	sp P80513	4
Nuclear migration protein nudC	sp Q17QG2	4
Lysosomal alpha-mannosidase	sp Q29451	4
Complement C2	sp Q3SYW2	4
Acyl-CoA-binding domain-containing protein 7	sp Q3SZF0	4
Mitochondrial fission 1 protein	sp Q3T0I5	4
Tetratricopeptide repeat protein 1	sp Q3ZBR5	4
Ran-binding protein 3	sp Q4R4T9	4
Protein DJ-1	sp Q5E946	4
Talin-1	sp Q9Y490	4
Vesicle-associated membrane protein-associated protein B	sp A2VDZ9	4
Decorin	sp P21793	4
Myotrophin	sp P58546	4
Transcription factor BTF3 homolog 4	sp Q2KIY7	4
Cysteine and glycine-rich protein 1	sp Q3MHY1	4
Probable D-tyrosyl-tRNA(Tyr) deacylase 2	sp Q96FN9	4
Dynein light chain 1, axonemal	sp Q2KID4	4
DnaJ homolog subfamily A member 2	sp O35824	4
Beta-2-glycoprotein 1	sp P33703	4
Dynein light chain 2, cytoplasmic	sp Q3MHR3	4
40S ribosomal protein S3a	sp B0KW94	3
Ubiquitin-like protein 4A	sp B2KIK3	3
Synaptotagmin-1	sp P21579	3
Sulfhydryl oxidase 1	sp O00391	3
Pigment epithelium-derived factor	sp Q95121	3
Prelamin-A/C	sp P02545	3
Large proline-rich protein BAT3	sp A5D9M6	3
Heterogeneous nuclear ribonucleoprotein Q	sp O60506	3
Chromobox protein homolog 3	sp Q13185	3
Protein S100-A13	sp P79342	3
Collagen alpha-1(XVIII) chain	sp P39061	3
Inter-alpha-trypsin inhibitor heavy chain H5	sp A2VE29	3
Agrin	sp O00468	3
Hemoglobin subunit alpha	sp P01966	3
Cofilin-1	sp P10668	3
Pleiotrophin	sp P21782	3
Catenin alpha-1	sp P26231	3
Annexin A11	sp P27214	3
Cholinesterase	sp P32749	3
Antithrombin-III	sp P41361	3
Arylsulfatase K	sp Q148F3	3
Poly(rC)-binding protein 2	sp Q15366	3

High mobility group protein B3	sp Q32L31	3
PDZ and LIM domain protein 7	sp Q3SX40	3
Leucine zipper transcription factor-like protein 1	sp Q3ZBL4	3
Transcobalamin-2	sp Q9XSC9	3
40S ribosomal protein S20	sp A1XQU9	3
Follitropin subunit beta	sp P04837	3
Thrombospondin-1	sp P07996	3
Microtubule-associated protein 4	sp P27546	3
Gamma-aminobutyric acid receptor-associated protein-like 1	sp P60518	3
Retinoic acid receptor responder protein 2	sp Q29RS5	3
PDZ and LIM domain protein 4	sp Q3T005	3
Serine/arginine repetitive matrix protein 1	sp Q5ZMJ9	3
Insulin-like growth factor-binding protein 7	sp Q61581	3
Aspartyl-tRNA synthetase, cytoplasmic	sp Q3SYZ4	2
Heterochromatin protein 1-binding protein 3	sp Q08DU9	2
Epsin-2	sp O95208	2
Insulin-like growth factor-binding protein 3	sp P20959	2
PC4 and SFRS1-interacting protein	sp Q8MJG1	2
THO complex subunit 4	sp B5FXN8	2
Vacuolar protein sorting-associated protein 4B	sp O75351	2
Ubiquitin-conjugating enzyme E2-17 kDa	sp P25867	2
Endoplasmic reticulum resident protein 29	sp P30040	2
Activated RNA polymerase II transcriptional coactivator p15	sp P53999	2
Tissue factor pathway inhibitor 2	sp Q7YRQ8	2
Glycogen phosphorylase, liver form	sp Q0VCM4	2
Protein S100-A10	sp P04163	2
Cystatin-B	sp P25417	2
Small nuclear ribonucleoprotein Sm D3	sp P62318	2
Ribosomal L1 domain-containing protein 1	sp A4FV97	2
Transforming growth factor-beta-induced protein ig-h3	sp O11780	2
NSFL1 cofactor p47	sp O35987	2
PDZ and LIM domain protein 3	sp O70209	2
Epsin-1	sp O88339	2
Kinectin	sp O97961	2
Histone H1.4 (Fragment)	sp P02252	2
Collagen alpha-2(I) chain	sp P02465	2
Fibronectin	sp P02751	2
Sodium/potassium-transport. ATPase subunit beta-1	sp P05028	2
Collagen alpha-3(VI) chain	sp P12111	2
ATP synthase subunit O, mitochondrial	sp P13621	2
Lamin-B1	sp P14731	2
Oxysterol-binding protein 1	sp P16258	2
Mimecan	sp P19879	2
Secretogranin-1	sp P23389	2

Tenascin	sp P24821	2
Peptidyl-prolyl cis-trans isomerase FKBP2	sp P26885	2
Elongation factor 1-beta	sp P34826	2
Protein S100-A4	sp P35466	2
60S ribosomal protein L30	sp P58372	2
DnaJ homolog subfamily B member 11	sp P81999	2
Eukaryotic translation initiation factor 3 subunit J	sp Q0VCU8	2
Ubiquitin-conjugating enzyme E2 variant 1	sp Q13404	2
ELAV-like protein 1	sp Q15717	2
CB1 cannabinoid receptor-interacting protein 1	sp Q17QM9	2
Heat shock 70 kDa protein 1B	sp Q27965	2
Collagen alpha-1(XII) chain (Fragment)	sp Q28902	2
Suppressor of G2 allele of SKP1 homolog	sp Q2KIK0	2
Vitamin D-binding protein	sp Q3MHN5	2
Hepatoma-derived growth factor-related protein 2	sp Q3UMU9	2
PDZ and LIM domain protein 5	sp Q96HC4	2
Cadherin EGF LAG seven-pass G-type receptor 2	sp Q9HCU4	2

## Additional file 3

Table 3 (complete, including accession number)

Complete list of identified proteins in the rhTSH preparation (Thyrogen, Genzyme GmbH; 2 lot numbers) by database search following mass spectrometry. Numbers of assigned spectra are given and the minimal sequence coverage was set to a minimum of 10%. Contaminations like keratin or trypsin, which were detected in bTSH as well as in rhTSH were excluded from analyses.

Identified Proteins	Accession Number	LOT A8035H40	LOT A8063H19
Thyrotropin subunit beta	sp P01222	116	126
Glycoprotein hormones alpha chain	sp P01215	25	26
Glyceraldehyde-3-phosphate dehydrogenase	sp P04406	6	1
Calmodulin-like protein 5	sp Q9NZT1	5	1
Caspase-14	sp P31944	4	1
Serpin B12	sp Q96P63	5	0
Galectin-7	sp P47929	2	0



## *Curriculum vitae*

### **Persönliche Daten**

Name, Vorname:	Schäfer Sandra
Anschrift:	Widenbachstrasse 24 8880 Walenstadt
Telefon:	079 / 613 14 50
E-Mail:	sandra_schaefer@gmx.ch
Geburtsdatum:	02.06.1983 in Niederbipp
Heimatort:	Lupsingen
Nationalität:	Schweiz
Zivilstand:	ledig

### **Ausbildung**

09.1998 - 09.2002	Gymnasium an der Kantonsschule in Solothurn Maturitätsprofil N (Mathematik und Naturwissenschaften)
09.2003 - 10.2009	Studium der Veterinärmedizin an der Vetsuisse Fakultät, Universität Zürich
10.2009	Diplom vet med an der Vetsuisse Fakultät, Universität Zürich
10.2009 - 03.2014	Doktorarbeit an der Vetsuisse Fakultät, Universität Zürich, Klinik für Kleintiermedizin
07.2011 - 07.2012	Internship an der Vetsuisse Fakultät, Universität Zürich, Klinik für Kleintiermedizin
07.2012 - heute	Residency ACVIM an der Vetsuisse Fakultät, Universität Zürich, Klinik für Kleintiermedizin

---