



**University of
Zurich^{UZH}**

**Zurich Open Repository and
Archive**

University of Zurich
University Library
Strickhofstrasse 39
CH-8057 Zurich
www.zora.uzh.ch

Year: 2018

All cats are gray in the dark: enrichment/depletion approaches for biomarker discovery on felis catus plasma

Carcangiu, L ; Pisanu, S ; Tore, S ; Addis, M F ; Zini, Eric ; Uzzau, S ; Pagnozzi, D

Abstract: In veterinary medicine, assay performance is often affected by the lack of species-specific diagnostic tools. Reliable biomarkers might be identified by investigating biological fluids of the species of interest, but protein sequence databases are often incomplete and human-specific devices for reducing sample complexity might fail when applied to animal plasma. Here, seven commercial methods based on different capturing agents (anti-human antibodies, affinity ligands, mixture of antibodies and ligands, and combinatorial peptide ligand libraries) are applied to cat plasma and evaluated in terms of yield, identified proteins/ peptides, and relative abundance by high-resolution shotgun proteomics and label-free quantitation. As a result, anti-human antibody-based methods are unsatisfactory. Most fail in reducing albumin and immunoglobulins, and some lead to a substantial removal of other highly abundant proteins, probably because of nonspecific interactions. A protein A/dye ligand-based method is efficient in reducing immunoglobulins, fibrinogen, and apolipoprotein A1 and A2, but not albumin, and protein identifications do not increase. Only peptide ligand libraries flatten the dynamic range, and increased protein identification (59.0%). Albumin and immunoglobulins are successfully depleted (60.7% and 35.9%, respectively). Although further studies will be required for reinforcing our observations, this work can provide a useful guide for cat plasma pretreatment in biomarker discovery studies.

DOI: <https://doi.org/10.1002/pmic.201800191>

Posted at the Zurich Open Repository and Archive, University of Zurich

ZORA URL: <https://doi.org/10.5167/uzh-159351>

Journal Article

Accepted Version

Originally published at:

Carcangiu, L; Pisanu, S; Tore, S; Addis, M F; Zini, Eric; Uzzau, S; Pagnozzi, D (2018). All cats are gray in the dark: enrichment/depletion approaches for biomarker discovery on felis catus plasma. *Proteomics*, 18(20):e1800191.

DOI: <https://doi.org/10.1002/pmic.201800191>

1 *Technical brief*

2

3 **All cats are grey in the dark: enrichment/depletion approaches for biomarker**
4 **discovery on *Felis catus* plasma.**

5

6 Running title: Enriching/depleting cat plasma samples

7

8 Laura Carcangiu¹, Salvatore Pisanu², Silvia Tore³, Maria Filippa Addis^{2,a}, Eric Zini^{1,4,5}, Sergio

9 Uzzau^{2,b}, Daniela Pagnozzi^{2*}

10

11 *1 Istituto Veterinario di Novara, Granozzo con Monticello, Novara, Italy*

12 *2 Porto Conte Ricerche, Science and Technology Park of Sardinia, Tramariglio, Alghero (Sassari),*

13 *Italy*

14 *3 Sardegna Ricerche, Piscina Manna, Pula, Cagliari, Italy*

15 *4 Clinic for Small Animal Internal Medicine, Vetsuisse Faculty, University of Zurich, Zurich,*

16 *Switzerland*

17 *5 Department of Animal Medicine, Production and Health, University of Padova, Legnaro, Italy.*

18 *a Present address: Department of Veterinary Medicine, University of Milan, Italy*

19 *b Permanent address: Department of Biomedical Sciences, University of Sassari, Sassari, Italy*

*Correspondence: Daniela Pagnozzi, PhD, Porto Conte Ricerche Srl, S.P. 55 Porto Conte - Capo

Caccia, km 8,400 Loc. Tramariglio - 07041 Alghero (SS) – Italy. Tel: (+39) 079 998 400 – Fax:

(+39) 079 998 567 (pagnozzi@portocontericerche.it).

20 Keywords: biomarker discovery; depletion; enrichment; plasma; shotgun proteomics.

21 Total number of words: 2406.

22 **Abstract**

23 In veterinary medicine, assay performance is often affected by the lack of species-specific
24 diagnostic tools. Reliable biomarkers might be identified by investigating biological fluids of the
25 species of interest, but protein sequence databases are often incomplete and human-specific devices
26 for reducing sample complexity might fail when applied to animal plasma. Here, seven commercial
27 methods, based on different capturing agents (anti-human antibodies, affinity ligands, mixture of
28 antibodies and ligands, and combinatorial peptide ligand libraries), were applied to cat plasma and
29 evaluated in terms of yield, identified proteins/ peptides, and relative abundance by high resolution
30 shotgun proteomics and label-free quantitation. As a result, anti-human antibody-based methods
31 were unsatisfactory. Most failed in reducing albumin and immunoglobulins, and some led to a
32 substantial removal of other highly abundant proteins, probably because of nonspecific interactions.
33 A protein A/dye ligand-based method was efficient in reducing immunoglobulins, fibrinogen and
34 apolipoprotein A1 and A2, but not albumin, and protein identifications did not increase. Only
35 peptide ligand libraries flattened the dynamic range, and increased protein identification (59.0%).
36 Albumin and immunoglobulins were successfully depleted (60.7% and 35.9%, respectively).
37 Although further studies will be required for reinforcing our observations, this work can provide a
38 useful guide for cat plasma pretreatment in biomarker discovery studies.

39

40 Farm animal and pet welfare is the main goal of veterinary medicine for both economic and
41 affection reasons. Nowadays, the performance of animal diagnostic serology is often sub-optimal,
42 mostly due to the low availability of commercial species-specific immunoassays [1]. Even though
43 clinical care may capitalize from what has already been discovered by human system biology,
44 dedicated biomarker discovery and validation efforts would be needed for improving veterinary
45 diagnostics [2]. Nevertheless, compared to humans, biomarker discovery studies directly performed
46 on animal specimens lag considerably behind [3-5]. This might be partly due to specific difficulties
47 in treating samples and to the limited availability of annotated protein sequence databases [6]. One
48 of the main problems in proteomic analysis of biological specimens such as plasma, serum, and
49 urine, is the presence of a huge dynamic molecular range of protein concentration [7-9]. The
50 availability of commercial products aimed at depleting highly abundant protein and enriching the
51 low abundant ones generally enables a more efficient and less challenging biomarker discovery
52 process. However, most commercial solutions based on immunoaffinity criteria have been
53 developed for human specimens, and only a limited amount of products are specific for, or
54 compatible with some animal species [10-13]. Other products are composed by a mixture of affinity
55 ligands (protein G or protein A, dye ligands, peptide libraries), theoretically allowing their use on a
56 wider range of organisms [14-16]. To our knowledge, no antibody-based depletion kit has been
57 declared by the manufacturers as suitable for feline samples, and only a limited number of papers
58 report proteomic studies on cat biofluids [17], none of them applying procedures for
59 depletion/enrichment of highly abundant proteins. Recently, a comparative evaluation of seven
60 commercial products was performed by our group on a human serum sample by means of a shotgun
61 proteomics approach [18, 19]. Here, we aimed at evaluating the same methods on a cat plasma
62 sample to gather information on their depletion performances and to provide a selection guide for
63 plasma biomarker discovery studies in this species.

64 Different capturing criteria were involved: i) immunoaffinity with antibodies directed against
65 specific (human) proteins, ii) protein G and antibodies directed against human serum albumin
66 (HSA), iii) protein A and cibacron blue dye-ligand, and iv) affinity for a combinatorial peptide
67 ligand library. All products were applied to a pool of plasma samples and the enriched proteins
68 were evaluated by SDS-PAGE, high resolution shotgun MS analysis and label-free quantitation.
69 Plasma samples from four cats, were homogeneously pooled and stored at -80°C until use. Pool
70 protein concentration was 63 mg/mL based on the BCA assay (Thermo Fisher Scientific, CA).
71 Depletions and enrichments were performed according to the product manuals with the sample
72 volumes recommended by the manufacturers [18]. Products were based on four different capturing
73 agents: antibodies (Qproteome Albumin/IgG Depletion kit, “Qproteome”, Qiagen, Sweden;
74 ProteoPrep Immunoaffinity Albumin and IgG Depletion Kit, “ProteoPrep”, Sigma-Aldrich, MO;
75 Top 2 Abundant Protein Depletion Spin Columns “Top 2”, and Top 12 Abundant Protein Depletion
76 Spin Columns “Top 12”, Thermo Fisher Scientific), specific ligands (Albumin/IgG Removal,
77 “CibB-A”, Thermo Fisher Scientific), mixture of antibodies and ligands (Albumin and IgG
78 Depletion SpinTrap, “SpinTrap”, GE Healthcare, Sweden), and combinatorial peptide ligand
79 libraries (ProteoMiner beads, “ProteoMiner”, Bio-Rad, CA), respectively. Table S1 summarizes the
80 selected products and their relative features. Two technical replicates were performed for each
81 procedure. A non-reducing Laemmli buffer was added to 10 µg of each depleted/enriched fraction
82 and SDS-PAGE was carried out on any kD precast polyacrylamide gels using the Mini-Protean
83 system (Bio-Rad). Filter aided sample preparation (FASP) was performed on all protein samples to
84 obtain tryptic peptide mixtures; LC-MS/MS analysis was carried out on a Q Exactive mass
85 spectrometer interfaced with an UltiMate 3000 RSLCnanoLC system (Thermo Fisher Scientific), by
86 performing a linear gradient of 245 min from 5% to 37.5% of eluent B (0.1% formic acid in 80%
87 acetonitrile) in eluent A (0.1% formic acid) at a flow rate of 250 nL/min. MS data were acquired
88 using a data-dependent top12 method, as previously described [18, 21]. Proteome Discoverer

89 (version 1.4; Thermo Fisher Scientific) was used for protein identification with Sequest-HT as
90 search engine. Raw files were processed with the following parameters: database *Felis silvestris*
91 *catus*, downloaded from UniProtKB_TrEMBL (release 2017 08); enzyme trypsin, with two missed
92 cleavages allowed; precursor mass tolerance 10 ppm; MS/MS tolerance 0.02 Da; charge states +2,
93 +3, and +4; cysteine carbamidomethylation as static modification and methionine oxidation as
94 dynamic modification. The percolator algorithm was used for protein significance and for peptide
95 validation (FDR<0.01%). Uncharacterized proteins were identified by Basic Local Alignment
96 Search Tool (BLAST) homology search.

97 Relative protein abundance was assessed by label-free quantitation based on the spectral counting
98 (SpC) approach, calculating the normalized spectral abundance factor (NSAF) according to Old et
99 al. [20]. All statistical analyses were performed using *t*-test (p-value ≤ 0.05).

100 According to the manufacturers suggestions, 10 μ L of plasma were processed by Top 2, Top 12 and
101 CibB-A columns, approximately corresponding to 0.63 mg of proteins; 25 μ L of sample were
102 incubated with SpinTrap, Qproteome and ProteoPrep columns, equivalent to 1.58 mg of proteins,
103 and a 40 μ L aliquot, accounting for 2.52 mg of proteins, was incubated with the appropriate volume
104 of ProteoMiner beads, respecting the volume/volume ratio suggested by the manufacturer for
105 serum/plasma samples. To evaluate protein yields for each product, the concentrations of
106 depleted/enriched fraction were estimated. Unexpectedly, total recovery was very high in most of
107 the applied protocols; more in detail, the protein content of recovered fractions spanned from a
108 minimum of 31% for Top 12, to 33% for CibB-A, 35% for ProteoPrep, 39% for Top 2, up to a
109 maximum of 48% for SpinTrap and Qproteome, when compared to the untreated plasma sample.
110 Only ProteoMiner led to a very low yield of 1.8%, suggesting a depletion of a large amount of
111 proteins. These yields were very different from the analogues calculated on human serum [18],
112 revealing that both immuno- and chemical-affinity based products, except for the combinatorial
113 peptide ligand library, have weaker interactions with cat plasma proteins.

114 To visualize the overall effect of each method, 10 μ g aliquots of untreated plasma, of the
115 depleted/enriched fractions and of the residual proteins were analyzed by SDS-PAGE. As shown in
116 Figure 1, as a result of all procedures but ProteoMiner, the protein profile (red) remained very
117 similar to the untreated sample (purple). However, a partial removal of highly abundant proteins did
118 occur, as demonstrated by the residual protein fractions collected during each procedure (Figure 1,
119 blue), which revealed a depletion of feline serum albumin (FSA, gel band at about 50 kDa) and
120 immunoglobulins (Igs, gel bands from the top of the gel to about 100 kDa), even though in different
121 amounts. After tryptic digestion, untreated plasma and depleted/enriched fractions were analyzed by
122 a long-gradient LC-MS/MS to assess protein numbers and identities, and the relative abundance
123 was estimated by label-free quantitation. Table 1 and Supplementary Table S2 summarize MS
124 results. The highest number of proteins and peptides were identified after ProteoMiner enrichment,
125 showing an increase of 59.0% and 24.9% when compared to untreated plasma, respectively. The
126 other procedures led to a number of identified proteins quite similar to the undepleted plasma. Only
127 CibB-A reduced the number of identified proteins and peptides, probably due to nonspecific
128 interactions with many plasma proteins.

129 Further, to assess the effect of treatments on most abundant proteins, their variation in comparison
130 to the untreated sample was calculated. Considering that the sequence identity between HSA and
131 FSA is 82.4%, it would be reasonable to expect an efficient removal of this protein with an anti-
132 HSA antibody, as reported by others for dog serum albumin (80.1% identity with HSA) [13].
133 Surprisingly, FSA removal was efficient only with the ProteoMiner approach, with a depletion of
134 60.7%; FSA content was slightly reduced by Top 2, whereas its amount was even increased by all
135 the other products. CibB-A, despite exploiting a dye ligand bait supposed not to be species-specific,
136 provided a 2.2-fold increase in the relative FSA content (Figure 2A). Regarding Igs, all the
137 identifications referable to Igs were considered and the overall variation was calculated (Figure 2B).
138 The best result was obtained by CibB-A, leading to a 59.7% of total Igs decrease, suggesting that

139 protein A has a higher affinity to cat Igs than protein G and anti-Igs antibodies. In general, all
140 procedures allowed only a modest removal of Igs. Interestingly, results were similar to those seen
141 with human serum [18] on other highly abundant proteins; more in detail, upon ProteoMiner
142 treatment, serotransferrin, α 2-macroglobulin, α 1-antitrypsin, haptoglobin, and α 1-acid glycoprotein
143 (Accession: M3WBQ5, M3WMA9, M3WCX1, M3W4S3, and Q6KCA0, respectively) were
144 decreased, whereas fibrinogen α -chain and fibrinogen β -chain (Accession: M3W022, and M3WII3)
145 were increased; likewise, CibB-A led to a considerable reduction of fibrinogen. In addition, CibB-A
146 induced a substantial decrease of apolipoprotein A1 and A2 (Accession: M3WPG6, and M3WN87)
147 and an increase of serotransferrin, α 2-macroglobulin, α 1-antitrypsin, haptoglobin, and α 1-acid
148 glycoprotein contents. The other products, apart from few variations consisting in a slight decrease
149 of fibrinogen α -chain and fibrinogen β -chain by Top 12, an increase of α 1-antitrypsin by both Top 2
150 and Top 12, and of α 1-acid glycoprotein by SpinTrap and QProteome, did not have an effect on the
151 most abundant proteins. Further, the effect on low abundant proteins was evaluated; a total of 302
152 proteins were identified only by ProteoMiner approach (Supplementary Table S3). Among them,
153 many human homologous proteins have been reported to be present in ng/mL (tenascin X, [22],
154 thrombomodulin [23] and others) to pg/mL (inhibin B, [24]) concentrations in serum samples from
155 normal subjects.

156 In conclusion, our data demonstrate that anti-human antibody-based methods are not suitable for cat
157 plasma samples; on the contrary, CibB-A was very efficient in reducing Igs, fibrinogen and
158 apolipoprotein A1 and A2; noteworthy, ProteoMiner allowed flattening of the sample dynamic
159 range and a significant increase in the number of identified proteins. It is likely that the binding
160 affinities of anti-human antibodies, protein G and cibacron blue dye-ligand are not adequate for cat
161 proteins removal, even in case of high sequence identity with the homologous human protein,
162 maybe due to a different protein domain exposure, or to the absence/presence of post-translational
163 modifications that impairs the interactions. Otherwise, the high heterogeneity of ligands provided

164 by ProteoMiner beads (several millions of hexapeptide sequences) might favor the chance of
165 establishing interactions . Moreover, it should be underlined that all procedures, but ProteoMiner
166 beads, provide for subsequent analysis a flow-through fraction, that may still be rich in target
167 proteins due to their overabundance compared to the limited availability of binding sites;
168 consequently, a fine optimization of the sample amount to be loaded should be performed before
169 treatments. Although further studies carried out on a higher number of samples will be required to
170 reinforce our observations, this work provides a useful guide for selecting depletion procedures in
171 biomarker discovery studies on cat plasma.

172

173 This work was supported by Regione Autonoma della Sardegna, with the grants “Diagnostici per
174 ANimali Da REddito e da Affezione”, D.AN.D.RE.A., and “art. 9 LR 20/2015”.

175 The authors declare no conflict of interest.

176 The authors thank Drs Tiziana Cubeddu and Giovanna Casalloni for providing plasma samples.

177

178 **References**

- 179 [1] T. Warnken, K Huber, K. Feige, BMC Vet. Res. 2016, 12, 196.
- 180 [2] R.E. Moore, J Kirwan, M.K. Doherty, P.D. Whitfield, Biomark. Insights 2007, 10 185.
- 181 [3] L. Galeandro, N.S. Sieber-Ruckstuhl, B. Riond, S. Hartnack, R. Hofmann-Lehmann, C.E.
- 182 Reusch, F.S. Boretta, J. Vet. Intern. Med. 2014, 28, 1433.
- 183 [4] R. Tanoue, I. Kume, Y. Yamamoto, K. Takaguchi, K. Nomiya, S. Tanabe, T. Kunisue, J.
- 184 Chromatogr. A 2018, 2, 30.
- 185 [5] A. Tanca, D Pagnozzi, G.P. Burrai, M. Polinas, S Uzzau, E Antuofermo, M.F. Addis, J.
- 186 Proteomics 2012, 21, 561.
- 187 [6] F. Di Girolamo, A. D'Amato, I. Lante, F. Signore, M. Muraca, L. Putignani, Int. J. Mol. Sci.
- 188 2014, 15, 15396.
- 189 [7] G. Biosa, M.F. Addis, A. Tanca, S. Pisanu, T. Roggio, S. Uzzau, D. Pagnozzi, J. Proteomics
- 190 2011, 75, 93.
- 191 [8] A.H. Zhang, H. Sun, G.L. Yan, Y. Han, X.J. Wang, Appl. Biochem. Biotechnol. 2013, 170, 774.
- 192 [9] S. Filip, K. Vougas, J. Zoidakis, A. Latosinska, W. Mullen, G. Spasovski, H. Mischak, A.
- 193 Vlahou, J. Jankowski, PLoS One 2015, 10, e0133773.
- 194 [10] R. Günther, E. Krause, M. Schümann, J. Ausseil, J.M. Heard, I.E. Blasig, R.F. Haseloff, Fluids
- 195 Barriers CNS 2014, 23, 14.
- 196 [11] K. de Morais-Zani, K.F. Grego, A.S. Tanaka, A.M. Tanaka-Azevedo, J. Biomol. Tech. 2011,
- 197 22, 67.
- 198 [12] D.R. Haudenschild, A. Eldridge, P.J. Lein, B.A. Chromy, Biochem. Biophys. Res. Commun.
- 199 2014, 455, 84.
- 200 [13] M. Sundberg, J. Bergquist, M. Ramström, Biochem. Biophys. Rep. 2015, 3, 68.
- 201 [14] K. Björhall, T. Miliotis, P. Davidsson, Proteomics 2005, 5, 307.

202 [15] K. Leecharoenkiat, W. Sornjai, K. Khungwanmaythawee, A. Paemanee, C. Chaichana, S.
203 Roytrakul, S. Fucharoen, S. Svasti, D.R. Smith, Dis. Markers 2014, 2014, 340214.

204 [16] A. D'Amato, A. Bachi, E. Fasoli, E. Boschetti, G. Peltre, H. Sénéchal, P.G. Righetti, J.
205 Proteome Res. 2009, 8, 3925.

206 [17] M.D. Meachem, E.R. Snead, B.A. Kidney, M.L. Jackson, R. Dickinson, V. Larson, E. Simko,
207 Can. J. Vet. Res. 2015, 79, 184.

208 [18] S. Pisanu, G. Biosa, L. Carcangiu, S. Uzzau, D. Pagnozzi, Talanta, 2018, 185, 213.

209 [19] S. Pisanu, G. Biosa, L. Carcangiu, S. Uzzau, D. Pagnozzi, Data Brief 2018, 19, 1765.

210 [20] W.M. Old, K. Meyer-Arendt, L. Aveline-Wolf, K.G. Pierce, A. Mendoza, J.R. Sevinsky, K.A.
211 Resing, N.G. Ahn, Mol. Cell. Proteomics 2005, 4, 1487.

212 [21] D. Pagnozzi, F. Tamarozzi, A.M. Roggio, V. Tedde, M.F. Addis, S. Pisanu, G. Masu, C.
213 Santucci, A. Vola, A. Casulli, G. Masala, E. Brunetti, S. Uzzau, Clin. Infect. Dis. 2018, 66, 1342.

214 [22] K. Yamada, A. Watanabe, H. Takeshita, K. Matsumoto, Clin. Chim. Acta 2016, 459, 94.

215 [23] F. Califano, T. Giovanniello, P. Pantone, E. Campana, C. Parlapiano, F. Alegiani, G.M.
216 Vincentelli, P. Turchetti, Eur. Rev. Med. Pharmacol. Sci. 2000, 4, 59.

217 [24] T.W. Kelsey, A. Miles, R.T. Mitchell, R.A. Anderson, W.H. Wallace, PLoS One 2016, 11,
218 e0153843.

219

220 **Figure legends**

221 **Figure 1. SDS-PAGE of plasma fractions with and without treatments.** Protein profiles of
222 untreated plasma (purple), depleted/enriched fractions (red) and residual proteins (blue) after Top 2,
223 Top 12, SpinTrap, Qproteome, ProteoPrep, CibB-A, ProteoMiner. **U:** unbound protein fractions. **E:**
224 proteins eluted from the resins. **M:** Precision Plus molecular weight markers (Bio-Rad).

225

226 **Figure 2 Abundance (NSAF) variations of serum albumin (panel A) and total**
227 **immunoglobulins (panel B) after plasma treatment.** Plotted values represent the average
228 NSAF \pm SD. Asterisks indicate statistically significant differences between untreated and
229 depleted/enriched plasma according to t-test, with *p value \leq 0.05 and ** p value \leq 0.01.

230

231 **Table 1.** Numbers of identified proteins, peptides, PSMs and Search Inputs in untreated and
 232 depleted/enriched plasma.
 233

	Proteins	CV (%)	Peptides	CV (%)	PSMs	CV (%)	Search Inputs	CV (%)
Plasma	390 ± 8	1.99	4245 ± 63	1.48	13112 ± 86	0.65	80323 ± 86	0.11
CibB-A	264 ± 0**	0.00	2214 ± 30**	1.34	10455 ± 71**	0.68	70498 ± 1242**	1.76
ProteoPrep	353 ± 13	3.61	3932 ± 52*	1.33	12608 ± 66*	0.53	81329 ± 1114	1.37
SpinTrap	385 ± 16	4.04	4047 ± 65	1.61	13377 ± 213	1.60	78574 ± 1314	1.67
Top 12	395 ± 7	1.79	4303 ± 24	0.56	13690 ± 74*	0.54	79571 ± 1225	1.54
Qproteome	396 ± 6	1.61	3872 ± 69*	1.77	13488 ± 601	4.46	77057 ± 2502	3.25
Top 2	398 ± 7	1.78	4346 ± 13	0.29	13602 ± 71*	0.52	78296 ± 44**	0.06
ProteoMiner	620 ± 2**	0.34	5303 ± 42**	0.80	15408 ± 65**	0.42	77151 ± 619*	0.80

234 Numbers represent the average from technical replicates ± SD.
 235 CV: coefficient of variation. Significant differences between untreated and depleted/enriched plasma, according to *t*-test, are
 236 indicated as * (p value ≤ 0.05), or ** (p value ≤ 0.01).
 237